General Catalogue

# Poulten & Graf

**SUPERIOR LABORATORY PRODUCTS**

**Catalogue Edition 2015/2016**

Copyright, Layout, Text, Photos  
Poulten & Graf GmbH,  
Poulten & Graf Ltd.,  
Photo Studio Eilers, Stadtprozelten

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Copyright, Layout, Text, Fotos  
Poulten & Graf GmbH,  
Poulten & Graf Ltd.,  
Fotostudio Eilers, Stadtprozelten

Quality is  
our  
tradition.

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## Poulten & Graf Ltd.

Peak Works 1 Alfreds Way  
Barking Essex IG11 0AS, UK

VAT-No. GB749862867

t +44 (0)20 8594 4256

f +44 (0)20 8594 8419

VOLAC@poulten-graf.com

www.poulten-graf.co.uk

### Bank Account/Bankkonto

Allied Irish Bank (AIB)  
10 Berkeley Square  
London W1J 6AA

#### GBP Payments

Bank Sort Code: 23-83-97  
Account No.: 45708002

#### Euro Payments

Bank Sort Code 23-92-85  
Account Number 20141076

IBAN No.

GB92 AIBK 2392 8520 1410 76

BIC Code: AIB KGB 2L

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Karl-Carstens-Str. 10  
D-97877 Wertheim

VAT-Nr. DE812843788  
Steuer-Nr. 800822274

t +49 (0)9342 9229 0

f +49 (0)9342 9229 80

sales@poulten-graf.com

www.poulten-graf.de

### Bank Account/Bankkonto

Volksbank Main-Tauber eG, Wertheim  
IBAN: DE 91 6739 0000 0000 0732 02  
BIC-Code: GENO DE 61 WTH

Sparkasse TauberfrankenWertheim  
IBAN: DE 07 6735 2565 0003 0125 98  
SWIFT-BIC: SOLA DE S1TBB

## About us

The Poulten & Graf Group was formed, ready for the new Millennium, on 31st December 1999 by the bringing together of two long established family companies, Walter Graf & Co., based in Wertheim, Germany, and John Poulten Limited, based in Barking, England. The Group comprises Poulten & Graf Holdings GmbH with its two manufacturing subsidiaries, Poulten & Graf Limited and Poulten & Graf GmbH, both 100% owned by the Poulten & Graf Holdings GmbH.

The objective of the new enterprise was to improve the production capabilities of both the former companies and to increase market share throughout Europe and beyond. The Chairmen of the holding company are John Poulten, Claire Poulten and Hans Graf, with Hans Graf as the Managing Director of Poulten & Graf Limited and Poulten & Graf GmbH. Following the merger, production was rationalised and the product range extended to improve the breadth of service available to our customers.

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Barking

Additionally, strong relationships with our European partners have been established to further extend the business to create a prosperous organisation with a significant share of the glassware and liquid handling markets. Over the years constant efforts have been made to improve customer service and administration and to build on the existing strengths and reputation of the founding companies as manufacturers and suppliers of high quality glassware and liquid handling products.

The directors and employees of the Group companies are also keen to preserve the family character of the constituent companies and to provide secure employment for our staff as well as quality products for our customers.

The best known of the product ranges offered by the companies are FORTUNA® and  .

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## Mission Statement

It is Poulten & Graf's mission to design, manufacture, calibrate and market a range of laboratory glassware, laboratory liquid handling equipment, syringes, bottles and other laboratory products so as to produce:

- a profitable return on capital invested and employed
- fair working conditions for its employees
- value for money and excellent service for its customers.

It is specifically Poulten & Graf's mission to offer the highest levels of guaranteed volumetric accuracy when calibrating volumetric glassware and laboratory liquid handling devices. Products and services offered to Poulten & Graf's customers will be of the highest possible standards of quality taking into account its customers' needs and aspirations. It is Poulten & Graf's aim to work in true partnership with both its customers and suppliers to form lasting relationships.

It is Poulten & Graf's objective

- to foster warm business and personal relationships between its German and UK operating companies and to work together to rationalise activities for the best joint performance and efficiency
- to put customers first at all times
- to provide long term and enjoyable employment to its staff
- to work for the preservation of the environment.

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Poulten & Graf  
GmbH  
ISO 9001:2008

## Confidence

When you purchase our products, you give us your trust. We appreciate this trust and will do our best to supply only high-quality perfect products.

## Quality and Accuracy

FORTUNA® stands for more than 90 years of quality. Certification according to ISO 9001 and DIN EN ISO 13485 for medical products proves this claim of high quality and gives the customer a warranty of quality products.

The  mark comes from the English 'VOLumetric ACcuracy'. 'VOLAC'-products are manufactured according to BS EN ISO 9001. Volumetric glassware of class A accuracy is individually tested and printed. For class A an individual number gives a 100 % warranty of accuracy as well as full traceability and a level of reliability unobtainable under any other brand.

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There you will find further information about the Poulten & Graf history, about our products and the latest company news including our webshop (see next page).

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0049 9342 92290

## Our Webshop is always open!

At <http://shop.poulten-graf.de> or <http://shop.poulten-graf.co.uk> you can easily order the whole product range of Poulten & Graf.

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Experience perfection  
through a combination  
of quality, precision  
and design!

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## 02. Liquid Handling · Dosiertechnik

**No matter whether it is manual or automatic, aqueous solution or aggressive acid - we have the right device for your application!**

Due to our long experience in the production of manual bottle top dispensers and automatic dosing systems, our products are easy to use, safe and precise. They are indispensable in everyday life in the laboratory.

We are proud to present our wide product range, which extends from bottle top dispensers and microliter pipettors to an electronic titrating device and an automatic dosing system. These cover most laboratory applications and a number of chemical and pharmaceutical manufacturing uses as well. With the widest product range on the laboratory market we offer real choice and expertise to our customers.

We can answer your dosing questions. Please contact us!

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	Acids Säuren		HF HF	Alkaline Sol. Laugen	Aqueous Sol. wässrige Lösungen	Solvents Lösungsmittel
	weak	strong				
	schwach	stark				
 <b>FORTUNA® POLYFIX® 101.070*</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred
 <b>FORTUNA® POLYFIX® 101.071*</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred
 <b>FORTUNA® POLYFIX® 101.072*</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred
 <b>FORTUNA® OPTIFIX® BASIC</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred
 <b>FORTUNA® OPTIFIX® SOLVENT</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred
 <b>FORTUNA UNIVERSAL® OPTIFIX® HF</b>	Preferred	Preferred	Preferred	Preferred	Preferred	Preferred
 <b>FORTUNA® OPTIFIX® SAFETY</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred
 <b>FORTUNA® OPTIFIX® SAFETY S</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred
 <b>OPTIMAT® Pumps Pipetten</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred
 <b>VIAAC® SMART Pipettor</b>	Preferred	Preferred	Not recommended	Preferred	Preferred	Preferred

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\* all media against which polypropylene is resistant.

For all other liquids please contact us!

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## FORTUNA® POLYFIX® Dispenser

- FORTUNA® POLYFIX® - with glass piston**
- FORTUNA® POLYFIX® - with PTFE-coated piston**
- FORTUNA® POLYFIX® - with amber glass cylinder and glass piston**

**Due to the fact that PP is suitable for a lot of chemicals, the POLYFIX® dispenser can be used for most liquids.**


- visible dispensing procedure
- easy-to-use volume setting system
- valve block made of polypropylene with FIX adapter
- according to EN ISO 8655-5

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according to EN ISO 8655-5

#### 101.070

- with glass piston
- suitable for all chemicals (incl. solvents), that do not attack polypropylene and don't crystallize

#### 101.071

- with coated PTFE-piston
- suitable for acids and alkaline solutions that do not attack polypropylene

#### 101.072

- with amber glass cylinder and glass piston
- suitable for all light sensitive liquids that do not attack polypropylene

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The bottle is not included in delivery (Matching bottles you find on p. 260 ff)

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comes with 1 PTFE intake tube, 2 PP adaptors (GL 40 + 45), packing unit: 1

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### Replacement Parts - for 101.070

#### 101.077

Glass piston, complete with finger guard and volume setting system

#### 101.075

Valve block unit, complete with dispensing cylinder and discharge tube

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### - for 101.071

#### 101.078

PTFE coated piston, complete with finger guard and volume setting system

#### 101.075

Valve block unit, complete with dispensing cylinder and discharge tube

#### 101.089

Stopper caps for discharge tubes, natural, pack 10

#### 101.093

Intake tubes made of PTFE

### - for 101.072

#### 101.077

Glass piston, complete with finger guard and volume setting system

#### 101.079

Valve block unit, complete with amber dispensing cylinder and discharge tube

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#### 101.093

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## Replacement Parts for POLYFIX® Dispenser

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## FORTUNA® OPTIFIX® Dispenser

### Advantages at a glance

- Best safety by use of high grade material
- Visible dosing procedure
- No contamination by metal parts, e.g. springs
- High accuracy, according to EN ISO 8655-5
- Premium quality
- BASIC and SOLVENT - wide volume range up to 500 ml

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## 02. Liquid Handling · Dosiertechnik

### Dispensing Piston

- proven technology
- made of glass, solid PTFE-encased, exception: SOLVENT
- virtually indestructible

### Precision Dispensing Cylinder

- made of glass, dosing procedure visible (except type HF, made of PTFE)
- coated with plastic - protection against breakage
  - graduation clear readable
  - easy to dismantle and to clean

### Valve Block

- 100 % tightness
- no metal springs - no contamination of the dosing liquid

### Finger Protection Cap

- ergonomic type
- timeless design
- with colour-code-system

### Volume Setting System

- fast and easy volume setting
- additional fine adjustment
- for right- and left hand use
- easy to read scale

### FIX Adaptor

- to align the dispenser on a bottle
- made of PP for BASIC + SOLVENT
- made of PTFE for SAFETY +HF

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according to EN ISO 8655-5

**101.080**

**for dispensing light acids,  
bases, aqueous solutions  
and crystallizing liquids**

- with solid PTFE-encased dispensing piston
- dispensing cylinder made of glass
- visible dosing procedure
- no metal spring, no contamination of the dosing liquid
- volume up to 500 ml

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
**101.080**  
**200 - 500 ml**

**comes with 1 PTFE intake tube, 2 PP adaptors (GL 40 + 45), packing unit: 1**

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**101.081**

**for dispensing solvents**

- with glass piston
- dispensing cylinder made of glass
- visible dosing procedure
- no metal spring, no contamination of the dosing liquid
- volume up to 500 ml

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
according to EN ISO 8655-5

**101.106**

**For all liquids that can be handled with the OPTIFIX dispensers and for hydrofluoric acid**

- made of PTFE
- with PTFE-encased dispensing piston and PTFE-cylinder
- FIX adaptor made of solid PTFE
- no metal spring, no contamination of the dosing liquid
- volume up to 50 ml

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comes with 1 PTFE intake tube, 2 PTFE adaptors (GL 40 + 45), packing unit: 1

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according to EN ISO 8655-5

**101.107**

**for dispensing strong acids and aggressive chemicals**

- PTFE-encased dispensing piston with sealing lip
- dispensing cylinder made of glass
- visible dosing procedure
- FIX adaptor made of solid PTFE
- no metal spring, no contamination of the dosing liquid
- volume up to 100 ml

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according to EN ISO 8655-5

#### 101.107 S

**for dispensing strong fuming acids and very abrasive liquids**

- with safety stopcock and additional safety valve in the ventilation tube
- PTFE-encased dispensing piston with sealing lip
- dispensing cylinder made of glass
- visible dosing procedure
- FIX adaptor made of solid PTFE
- no metal spring, no contamination of the dosing liquid
- volume up to 100 ml

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### - for 101.080

#### 101.000

Piston BASIC, solid PTFE encased, with finger guard and volume setting system

#### 101.004

Valve block unit with dispensing cylinder

### - for 101.081

#### 101.001

Piston SOLVENT, made of glass, with finger guard and volume setting system

#### 101.004

Valve block unit with dispensing cylinder

### - for 101.107

#### 101.003

Piston SAFETY, solid PTFE encased, with sealing lip, finger guard and volume setting system

#### 101.004

Valve block unit with dispensing cylinder

### - for 101.107-..S

#### 101.003

Piston SAFETY S, solid PTFE encased, with sealing lip, finger guard and volume setting system

#### 101.004-S

Valve block unit with dispensing cylinder, with safety valve

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## Replacement Parts for OPTIFIX® Dispenser

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### Volume setting system

#### 101.000-KS

Volume setting system for BASIC, SOLVENT, HF, SAFETY, SAFETY S, 2-100 ml

#### 101.000-KM

Volume setting system for BASIC, SOLVENT, 200 - 500 ml

### Intake tubes

#### 101 093-36

Intake tubes for HF, 5 ml, made of PTFE

#### 101.093-38/-48/-51/-58

Intake tubes for BASIC, SOLVENT, HF (10+30 ml), SAFETY, SAFETY S, made of PTFE

#### 10.1433

Tubing, made of PTFE, sold in varying lengths (state length with order)

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### Discharge tube sets

#### 101.089

Stopper caps for discharge tubes, natural

#### 101.094

Discharge tube set (tube with stabiliser, fixing screw, cap)

#### 101.095

Discharge tube set with Safety-Stopcock (tube with stabiliser, 2 fixing screws, PTFE-stopcock with dispensing tip), completely interchangeable; prevents dripping after the dispensing procedure

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made of PTFE

#### 101.105

Humidity absorber, made of PTFE, for hygroscopic liquids, single pack

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**101.104**

for OPTIFIX® SAFETY,  
SAFETY S, OPTIMAT®

- flexible remote dispensing
- suitable for easy serial dispensing
- tubing and dispensing valve made of PTFE
- with new metal spring, chemical resistant
- prevents liquid from coming out and stops liquid flowing back

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**The WANDERER  
for flexible  
remote dispensing**

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### FORTUNA® OPTIFIX® titrier

#### Benefits

- Easy to use
- easy to read LCD-display
- visible dispensing process
- left or right hand operation
- no metal spring
- high precision due to high tech micro electronic system
- can be used to dispense aggressive liquids

#### Application

- in chemical laboratories
- in environmental laboratories
- in food laboratories
- for routine daily work

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## Technical Features

- with certification of conformity
- GS = safety approved
- cylinder cap. 50 ml, display max. 99.99 ml
- subdivision 0.01 ml, Accuracy =  $\pm 0.6\%$ , Tolerance =  $\pm 0.2\%$
- according to EN ISO 8655-5

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### LCD-display and buttons

- the size of the LCD-display figures is 10 mm
- the arrow direction indicators show clearly the FILL or TITration mode
- LCD-display and buttons are at 10° angle to ensure easy reading and operation
- the mode is switched automatically from FILL to TITration

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## The dispensing unit

- the high grade registered glass cylinder is amber stained to allow use with light sensitive materials
- dispensing process is always visible and, if necessary, air bubbles can be easily removed
- the precision dispensing piston, made of glass, is coated with solid PTFE
- excellent resistance both to aggressive acids and alkaline solutions

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## The valve block system

- made of pure PTFE
- can be directly mounted on bottles with GL 32
- has an additional ball valve in the ventilation tube

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according to EN ISO 8655-5

**101.102**

for dispensing strong acids,  
alkalis and aggressive liquids

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**101.102EEZ**

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**101.094-48T**

Discharge tube ‚titrier‘, with  
extra thin tip for titrating drop  
by drop

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### **FORTUNA® OPTIMAT® dispensing pump**

**101.126**

**for dispensing strong fuming acids and very abrasive liquids, single pack (intake and discharge tube sets are not included)**

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**Pumps for  
OPTIMAT®  
programmable dispenser:**

Packing unit: 1

**Pumpen für  
OPTIMAT®  
Dosierstation:**

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**101.126**

Volume Volumen	Grad. Grad.	A* <sub>≤ ±</sub> R* <sub>≤ ±</sub>	CV* <sub>≤ ±</sub> V* <sub>≤ ±</sub>	max.V max.V	Art.-No. Art.-Nr.
0 - 1 ml	0.05 ml	± 0.6 %	± 0.2 %	0,55 ml/s	101.126-20
0 - 2 ml	0.1 ml	± 0.6 %	± 0.2 %	1,90 ml/s	101.126-27
0 - 10 ml	0.2 ml	± 0.6 %	± 0.2 %	3,35 ml/s	101.126-38
0 - 30 ml	0.5 ml	± 0.6 %	± 0.2 %	9,60 ml/s	101.126-44
0 - 50 ml	1.0 ml	± 0.6 %	± 0.2 %	15,20 ml/s	101.126-48
0-100 ml	2.0 ml	± 0.6 %	± 0.2 %	24,30 ml/s	101.126-51
0-200 ml <sup>(1)</sup>	5.0 ml	± 0.6 %	± 0.2 %	31,30 ml/s	101.126-53
0-300 ml <sup>(1)</sup>	5.0 ml	± 0.6 %	± 0.2 %	31,80 ml/s	101.126-56

\* A=Accuracy, CV=Coefficient of variation \*R=Richtigkeit, V= Variationskoeffizient

<sup>(1)</sup> additional to ISO 8655 · in Ergänzung zur ISO 8655

### **Accessories for OPTIMAT® and OPTIMAT® dispensing pumps**

**101.132**  
Intake tube  
length 1 m

**101.134**  
Discharge tube set  
made of PTFE, length 0.75  
or 1.50 m

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**155.000-40**  
Serielles Verbindungskabel

**Tubes for  
150.000-3 OPTIMAT®  
programmable dispenser:**

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**Delivering PRECISION  
Delivering EXCELLENCE**

#### **The advantages at a glance:**

- High performance, light weight and easy to operate
- ergonomic design
- Fatigue free working with convenient finger rest
- Suitable for both left and right hand users with relaxed grip and a good balance
- Colour coding for easy identification
- Specially engineered plastics resulting in better protection against chemical and physical corrosion
- Fully autoclavable without disassembly
- Enhanced UV resistance
- Slim design allows pipette to operate even with small narrow tubes used in laboratory
- Supplied with individual QC and calibration report in accordance with ISO 8655
- CE IVD compliant
- 3 years warranty

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## Technical Details

- PTFE sealing is used to create precise air tightness to provide long term consistent performance
- Pipette handle has coating of Thermoplastic Elastomer which prevents transfer of body heat to the pipette - ensuring high accuracy even during continuous use
- Large and clear 4 digit-display giving smaller increments for wider volume selection options
- Volume display is always visible while pipetting
- Resistant free click stop counter
- Soft touch tip ejection
- New tip cone design for smooth fitting with most manufacturers' tips
- Easily accessible recalibration mechanism
- No chance of accidental change in calibration
- In house cleaning, repair and calibration possible without any special tools

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\* A=Accuracy, CV=Coefficient of variation

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**Packing unit: 1**

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**volac** SMART  
Micropipettor,  
8 Channels

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**Packing unit: 1**

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## 02. Liquid Handling · Dosiertechnik

**volac** SMART  
Micropipettor,  
12 Channels

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**Packing unit: 1**

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## 02. Liquid Handling · Dosiertechnik

### **volac** ECO Automatic Pipettes - from 0.1 ml to 10 ml

**E 1000**

#### The advantages at a glance:

- replace glass pipettes and pipette fillers
- remarkably accurate (to  $\pm 1.75\%$ )
- micrometer style volume setting by turning thumb press
- fully autoclavable at  $121^\circ\text{C}$
- easily dismantled, cleaned and serviced

### **volac** MiniPipet, fixed volume

**R570**

#### The advantages at a glance:

- low cost
- accuracies  $\pm 0.5\ \mu\text{l}$
- can be used with standard yellow  $200\ \mu\text{l}$  tips
- double action thumb press gives "blow out" or permits "reverse" pipetting
- fully autoclavable to  $121^\circ\text{C}$
- Individualised packing cards can be easily printed. Minimum orders 100 pieces (mixed volumes possible)

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**Packing unit:** 1

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## 02. Liquid Handling · Dosiertechnik

### A new range of calibrated Multifit micropipettor tips

#### The advantages at a glance:

- very good volumetric accuracy
- exceptional moulded quality
- excellent fit for accurate working
- suitable for SMART, ECO and other leading brands

### Disposable Micro-Pipettor Tips, calibrated

#### Index for type of packing:

- a. bulk pack
- b. open rack
- c. racked/boxed/  
not sterilised
- d. racked/boxed/sterilised

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**110.856-07**

#### The advantages at a glance:

- can be used with all glass or plastic pipettes between 0.1 – 200 ml
- easy single handed use due to light weight and ergonomic design
- easy change from fill to dispense via circuit feedback
- two fill and dispense speeds available
- long lasting batteries
- comes with hydrophobic PTFE membrane filter, 0.45 µm and 0.2 µm

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## 110.856-07

The  pipette filler  
(Art.-No. 110.856-07)  
is supplied complete with:

- charging base
- 1 x 0.45  $\mu\text{m}$  and 1 x 0.2  $\mu\text{m}$  hydrophobic PTFE-filters
- bench stand
- wall stand
- instruction manual

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**Packing unit:** 1

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**With the desalination unit LAB-ioN L 2 demineralised water can be produced in the simplest way.**

### **Advantages at a glance:**

- output approx. 30 - 40 l/h
- used cartridges are simply replaced
- connection tube to a 3/4" tap

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**Desalination unit LAB-IoN L 2 comes with:**

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**Packing unit: 1**

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**The Water Distillator produces distilled water from tap water. The mono distillate produced is free of germs and dissolved salts. It complies with DAB 8 regulations.**

**106.900**

Compact design, an automatic level regulation system, easy cleaning and a temperature sensor that switches off the distilling system in case of trouble, all argue for the Water Distillator.

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Precise and reliable –  
we know what is  
important to you.

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### **Poulten & Graf Volumetric Glassware - tested precision-made glassware**

Poulten & Graf offers a wide range of Volumetric Glassware under the trade marks FORTUNA® and  .

The whole range of Volumetric Glassware made by Poulten & Graf is produced according to British BS, European EN or international ISO standards.

There is also a specialised range made according to US Pharmacopoeia (USP) standards. In addition  cylinders and flasks are made to ASTM standards.

Poulten & Graf Volumetric Glassware is produced according to the highest quality standards. Control of raw materials, production and quality, work hand in hand and guarantee the high quality of the final product.

Poulten & Graf Volumetric Glassware, conformity certified, class A, is produced with a batch number and can be delivered with a batch certificate (except for Original FORTUNA® graduated and bulb pipettes with suction piston are delivered without certificate).

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## Graduation

We distinguish between two different graduation types:

- for Class A or AS products: complete ring marks
- for Class B products: short graduations

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## Calibration

Generally we differentiate between volumetric glassware which is calibrated to contain (IN) or calibrated to deliver (EX). Flasks, measuring and mixing cylinders are calibrated to contain (IN) so the liquid poured out corresponds with the volume indicated. Pipettes and burettes, however, are calibrated to deliver (EX). Provided the required waiting time is observed, the delivered liquid corresponds with the volume indication.

FORTUNA® and  glassware is individually calibrated (no mass production technique and no approximations).

**Special Note:** When users wish to check the accuracy of volumetric glassware, care should be taken not only to avoid parallax error, but also to avoid a false view of the meniscus. Use of a black ring or a clip approx. 1 mm below the mark will give a black reflection on the bottom of the meniscus and will dramatically improve the accuracy of the reading (see ISO 4787).

We can supply suitable black clips on request.

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### Colour marking

The colour markings on FORTUNA® and **volac** volumetric glassware are applied by silk screen printing. We only use diffusion or enamel colours. The diffusion colour - using amber stain - is completely absorbed into the glass surface during the firing process. It cannot be removed except by processes which remove the glass surface itself. It is therefore resistant to cleansing and even to very aggressive chemicals.

All other colours - e.g. blue, white, yellow, green - are enamel or melting colours, which are melted onto the glass surface at high temperature.

**volac** glassware is individually printed - both Class B and Class A. Particularly for graduated glass (graduated pipettes, burettes, cylinders) individual printing is the only way to ensure that the top and bottom of the scale matches individual calibration exactly.

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## Classes

In accordance with the ISO standard: FORTUNA®

Class	Certificate	Individually Numbered	Products
A	batch certificate	on request	cylinders, flasks
AS	batch certificate	on request	burettes, graduated and bulb pipettes

In accordance with the ASTM, ISO and USP standards: 

Class	Certificate	Individually Numbered	Products
B	no	no	burettes, graduated and bulb pipettes, cylinders
A	no	with or without individual certificate	burettes, graduated and bulb pipettes, cylinders, flasks
ASTM	single certificate	yes	cylinders, flasks
USP	single certificate	yes	burettes, graduated and bulb pipettes, cylinders, flasks

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## Why use volumetric glassware according to ASTM- or USP- standard?

In the United States and in Arabic Countries, as well as in many Asian Countries there is strong demand for volumetric glassware according to ASTM- or USP-standards.


Essentially, the tolerances of measuring devices made to these standards differ to a great extent from the European norm, ISO, and are much tighter.

If a company produces goods that have to comply with the American FDA import controls or come within their responsibility, they have to meet certain regulations (GMP).

GMP stands for Good Manufacturing Practice. This is a set of regulations, which ensures good quality, good production processes and appropriate production environments in the production of pharmaceutical products and chemical agents, as well as cosmetics, food and animal feed.

Amongst other things the focus lies in a complete and closed traceability of the products. In addition, the traceability and control of the manufacturing process (tracking), plays an important role.

To ensure the tracking of the volumetric glassware which Poulten & Graf produces, individual serial numbers and quality certificates are provided. Manufacturers, who do not follow the GMP and cannot ensure complete quality control, face warnings and penalties. As a last resort, producers, who aim to sell their products mainly in territory controlled by the FDA, can lose their permission to sell into the US market. To ensure maximum reliability, precision and traceability, choose **volac** volumetric glassware, produced according to ASTM- or USP- standard, with individual serial numbers and certificates.

### **volac** Volumetric Glassware Specifications

Volumetric Flasks Class A	ASTM and USP is tighter than ISO
Bulb Pipettes Class A	some USP are the same as ISO but some are tighter
Graduated Pipettes Class A	again some USP specifications are tighter than ISO but some are the same
Measuring Cylinders Class A	ASTM and USP specifications are unique
Burettes Class A	USP and ISO are identical

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### FORTUNA® Burettes

FORTUNA® burettes and automatic burettes are produced to the highest quality standards and only according to **Class AS (EN ISO 385)**.

**Reading:** With ring marks and Schellbach stripes, precise working is possible, reading the volume at the touch point of the two tips. Amber graduations are acid and alkali resistant. FORTUNA® burettes are calibrated to deliver ('EX'). Please remember the waiting time of 30 sec.

**Stopcocks:** FORTUNA® burettes are offered with straight stopcocks and glass or PTFE keys. FORTUNA® Automatic Burette, Pellet pattern are offered with lateral stopcock and glass or PTFE keys.

**Precision Burette Tips:** During the production of burettes we take great care with the burette tips so that the user can always titrate with the same size of drop.

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<sup>(3)</sup> class AS, without conformity approval

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### **FORTUNA® Burette, with straight stopcock and glass key**

#### **Class AS**

conformity approved,  
borosilicate glass 3.3

#### **1.008-...-02G, -04G**

Schellbach,  
blue or amber graduation

#### **1.008-...-06G**

glass amber stain diffused,  
white graduation

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### **FORTUNA® Burette, with straight stopcock and PTFE key**

**Class AS,**  
conformity approved,  
borosilicate glass 3.3

#### **1.008-..-02P, -04P**

Schellbach,  
blue or amber graduation

#### **1.008-..-06P**

glass amber stain diffused,  
white graduation

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
### **FORTUNA® Burette, Bang pattern,**

#### **Class AS,**

conformity approved (10 ml),  
borosilicate glass 3.3

#### **1.030**

with straight glass stopcock and glass  
key in upper stopcock, Schellbach,  
blue graduation

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**Packing unit: 1**

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### FORTUNA® Automatic Burette, Pellet pattern

**Class AS,**  
without conformity approval,  
automatic zero adjustment,  
borosilicate glass 3.3

#### 1.055..-02

without intermediate stopcock, with  
lateral PTFE stopcock, Schellbach,  
blue graduation

#### 1.055..-06

without intermediate stopcock, with  
lateral PTFE stopcock, amber stain  
diffused, white graduation

**Class AS,**  
conformity approved,  
automatic zero adjustment,  
borosilicate glass 3.3

#### 1.057..-02

with intermediate glass stopcock  
and lateral glass stopcock, Schell-  
bach, blue graduation

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
### FORTUNA® Automatic Burette, Pellet pattern

**Class AS,**  
without conformity approval,  
automatic zero adjustment,  
borosilicate glass 3.3

#### 1.059..-02

with intermediate PTFE stopcock and  
lateral PTFE stopcock, Schellbach,  
blue graduation

#### 1.059..-06

with intermediate PTFE stopcock and  
lateral PTFE stopcock, amber stain  
diffused, white graduation

### Reservoir bottle for automatic burettes

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
**Packing unit:** single pack,  
without reservoir

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### **FORTUNA® Automatic Burette, Dr. Schilling pattern**

#### **Class B**

automatic zero  
adjustment,  
soda lime glass

#### **1.068**

Schellbach,  
blue graduation

#### **1.069**

glass amber stain diffused,  
white graduation

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
#### **1.068**

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#### **1.069**

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comes mounted on a  
polyethylene bottle,  
with a plastic foot

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### **volac** Burettes


- manufactured in a range of sizes from 1 ml to 100 ml capacity
- **volac** burettes come in single packs
- produced according to BS ISO 385

### **volac** Burettes, Class A

- one cane form available: Class A standard
- featuring bold figures and fine graduation lines in white
- all Class A burettes are made of borosilicate glass
- have a flange at the top for ease of filling
- Class A burettes are available uncertified (suffix NC) or with a Works Accuracy Certificate (suffix WAC)
- all Class A burettes are individually calibrated and tested at 3 points (using registered tubing), all individually numbered and fully traceable
- removeable glass jets are inscribed with the number of the instrument for traceability
- all show their flow time inscribed on the instrument
- in addition Works Accuracy Certificate burettes come with a certificate specifying accuracy at three points

### **volac** Burettes, Class B

- three cane forms available: Class B standard, Class B Schellbach and Class B dispensing burettes
- each featuring bold figures and fine graduation lines in permanent amber
- Class B burettes are made of borosilicate glass or soda lime glass

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RB9803

**volac** Burette, with interchangeable glass key in glass barrel

### Class A

white graduation,  
borosilicate glass

**RB9803/NC**

without certificate

**RB9803/WAC**

with certificate

**volac** Burette, with interchangeable PTFE key in glass barrel

### Class A

white graduation,  
borosilicate glass

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
RB9807

**volac** Burette,  
2-way flow with  
interchangeable  
PTFE key in glass  
barrel

### Class A

white graduation,  
borosilicate glass

**RB9807/NC**

without certificate

**RB9807/WAC**

with certificate

**volac** Burette,  
with removeable  
PTFE-stopcock,  
with PTFE key,  
with glass tip

### Class A

white graduation,  
borosilicate glass

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
**volac** **Burette**  
**Class USP, A**

white graduation,  
borosilicate glass

**with interchangeable  
glass key in glass barrel**

**US9803/NC**

without certificate

**US9803/WAC**

with certificate

**with interchangeable  
PTFE key in glass barrel**

**US9806/NC**

without certificate

**US9806/WAC**

with certificate

**with 2-way flow, with  
interchangeable PTFE key  
in glass barrel**

**US9807/NC**

without certificate

**US9807/WAC**

with certificate

**with removable PTFE-  
stopcock, with PTFE key,  
with glass tip**

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**volac** **Burette, with interchangeable glass key in glass barrel**

**Class B**  
amber graduation,  
borosilicate glass

**R9103**

**Class B**  
Schellbach,  
borosilicate glass

**R9503**

**volac** **Burette, with interchangeable PTFE key in glass barrel**

**Class B**  
amber graduation,  
borosilicate glass

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**volac** Burette,  
2-way flow with  
interchangeable  
PTFE key in glass  
barrel

**Class B**

amber graduation,  
borosilicate glass

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braun graduiert,  
Borosilikatglas

**R9107**

**volac** Burette,  
2-way flow with  
interchangeable  
PTFE key in glass  
barrel, Schellbach

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**volac** Burette,  
with removeable  
PTFE-stopcock,  
with PTFE key,  
with glass tip

**Class B**  
amber graduation,  
borosilicate glass

**R9109**

**volac** Burette,  
with removeable  
PTFE-stopcock,  
with PTFE key,  
with glass tip

**Class B**  
amber graduation,  
borosilicate glass

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R200

**volac** **Dispensing Burette, with glass key**

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
**volac** Replacement  
complete PTFE stopcock

interchangeable complete  
PTFE stopcock,  
straight flow

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**R909**

**volac** Replacement  
glass burette jet,  
for removeable PP-  
stopcock

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**Packing unit: 1**

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### FORTUNA® Graduated Pipettes

FORTUNA® graduated pipettes are produced to the highest quality standards and are produced in Class AS. **(EN ISO 385)**.

**Reading:** FORTUNA® graduated pipettes are produced with ring marks (**Class AS**). Reading of the volume at the lowest point of the meniscus.

#### Technical Details:

- graduated pipettes are calibrated to deliver
- Colour-Code-System (ISO1769)
- short practical delivery time, openings uniformly calibrated
- tips and suction tubes reinforced
- fire polished facettes, bevels surface tempered
- amber graduations acid and alkali resistant
- indication of production batch (for conformity approved pipettes)
- soda lime glass

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\* additional to standard, without conformity approval

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\* additional to standard · ergänzend zur Norm

## FORTUNA® Graduated Pipette

total delivery, **Class B**,  
clear glass, with short  
graduation marks

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### FORTUNA® Graduated Pipette

total delivery, **Class AS**,  
conformity approved,  
clear glass, with ring marks

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### Original FORTUNA® Graduated Pipette with suction piston

clear glass, clear calibrated  
cylinders, without certificate

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## Graduated Pipettes

All  graduated pipettes offer exceptional standards of accuracy. Full ranges are available in **class B** and in **class A** (either uncertified or with a works certificate of accuracy).

- all graduation marks are in permanent amber stain
- colour codes are fused ceramic
- sizes 5 ml and above are constricted for cotton plugging
- all class A pipettes are individually tested at 5 places, numbered, timed and are fully traceable
- all  graduated pipettes satisfy both British and International Standards (**BS ISO 835**)
- all NC (uncertified) or WAC (certified) pipettes have a serial number and the flow time marked on them
- special Class A pipettes to USP are also available

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R377

### Graduated Pipette

**Class B**, amber stain graduation, individually calibrated and printed

#### R377 - Type 1

zero at top, capacity at shoulder, partial delivery

#### R376 - Type 2

zero at outlet, capacity at top, total delivery

#### R400 - Type 3

zero at top, capacity at outlet, total delivery, non blowout

#### R402 - Type 4

zero at top, capacity at outlet, total delivery, blowout

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**volac**

### Tissue Culture Pipette, short form

non British Standard, amber stain graduation

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R385

### Graduated Pipette

**Class A**, amber stain graduation, individually numbered

#### R386/NC - Type 1

zero at top, capacity at shoulder, partial delivery, without certificate

#### R386/WAC - Type 1

zero at top, capacity at shoulder, partial delivery, with certificate

#### R385/NC - Type 2

zero at outlet, capacity at top, total delivery, without certificate

#### R385/WAC - Type 2

zero at outlet, capacity at top, total delivery, with certificate

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**Packing unit: 5**

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### Graduated Pipette

**USP Class A**, amber stain graduation. Each pipette is individually numbered and fully traceable satisfying USP/ISO/BS/DIN accuracies.

#### US385 - Type 2

zero at outlet, capacity at top, total delivery, without certificate

#### US385/WAC - Type 2

zero at outlet, capacity at top, total delivery, with certificate

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### **FORTUNA® Pipette filler**

Is the most common pipetting aid for all types of pipettes. The filler has three glass ball valves which can be controlled by finger pressure.

**Functions:** releasing the air, filling the pipette, delivery of the liquid

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#### **1.170-04**

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### **PI-PUMP Pipette Filler**

easy-to-use design, with elastic step tip and drain valve

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
Packing unit: 1

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### FORTUNA® Bulb Pipettes

FORTUNA® bulb pipettes are produced to the highest quality standards. They are offered in **class AS** and **class AS conformity approved** ('A' stands for the highest quality specifications, 'S' for swift delivery). FORTUNA® bulb pipettes are calibrated to deliver ('EX') and are supplied with 1 or 2 ring marks. **(EN ISO 648)**.

**Reading:** Reading of the volume at the lowest point of the meniscus.

#### Technical Details:

- Colour-Code-System (ISO 1769)
- short practical delivery time, openings uniformly calibrated
- tips and suction tubes reinforced
- fire polished facettes, bevels surface tempered
- amber graduations acid and alkali resistant
- indication of production batch
- soda lime glass

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### 1.230 / 1.232 Bulb Pipettes, Class AS · Vollpipetten, Klasse AS

Capacity Inhalt	Tol. Tol.	Flow time Auslaufzeit	Total length Länge	Weight Gewicht
ml	± ml	s	mm	g
0.5	0.005	6 - 10	300	11
1	0.008	7 - 11	300	11
1.5	0.01	7 - 11	320	13
2	0.01	7 - 11	335	14
2.5	0.01	7 - 11	350	16
3	0.01	7 - 11	350	16
4	0.015	7 - 11	385	22
5	0.015	9 - 13	385	22
6	0.015	9 - 13	420	22
7	0.015	9 - 13	425	22
8	0.02	9 - 13	435	22
9	0.02	9 - 13	445	23
10	0.02	11 - 15	450	24
11	0.02	11 - 15	460	25
12	0.02	11 - 15	465	25
15	0.03	11 - 15	510	31
20	0.03	12 - 16	520	34
25	0.03	15 - 20	530	38
30	0.03	15 - 20	525	41
40	0.05	20 - 25	535	44
50	0.05	20 - 25	550	46
100	0.08	25 - 30	590	69

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### FORTUNA® Bulb Pipette

**Class AS**, with 1 mark,  
clear glass, conformity approved,  
with batch certificate

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### FORTUNA® Bulb Pipette

**Class AS**, with 2 marks,  
clear glass, conformity approved,  
with batch certificate

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### Original FORTUNA® Bulb Pipette with suction piston

clear glass, clear calibrated  
cylinders, without certificate

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### **volac** Bulb Pipettes

**volac** class B bulb pipettes are made to the same high standards and with the same one piece process as class A pipettes. The same amber stain markings, colour codes and superb accuracy of the jet are standard. All tips are 'Tufftip' treated to resist damage, all pipettes are calibrated to deliver.

**volac** class A bulb pipettes are guaranteed to give both volumes and flow times within class A limits. Suffix NC denotes verification only, suffix WAC denotes pipettes with a Works Accuracy Certificate.

All **volac** class A pipettes, whether NC or WAC, are individually calibrated and tested, individually numbered and fully traceable and each is individually marked with the time of flow to deliver.

All **volac** pipettes satisfy the International Standards (BS ISO 648) for accuracy and flow rate.

**volac** bulb pipettes according to USP Class A accuracy specifications are also available.

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**volac**

### Bulb Pipette

**Class B**, with 1 mark,  
amber stain graduation,  
individually calibrated

**R370**

**Class A**, with 1 mark,  
amber stain graduation,  
individually numbered

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### Bulb Pipette

**USP Class A**, with 1 mark, amber stain graduation, each piece is individually numbered and fully traceable satisfying USP/ISO/DIN/BS accuracies.

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### FORTUNA® Cylinders

FORTUNA® graduated and mixing cylinders (borosilicate glass 3.3) are 'Made in Germany' under the highest quality standards. They are produced in **Class B** or **Class A** (conformity approved). (**EN ISO 4788**).

**Reading:** FORTUNA® graduated and mixing cylinders are calibrated to contain ('IN'), with short graduation marks (Class B) or ring marks (Class A-DIN). Reading of the volume is at the lowest point of the meniscus.

#### Technical Details:

- delivered with hexagonal base or with plastic feet (autoclavable up to 130 °C),
- with tolerance indications
- amber graduations acid and alkali resistant
- indication of production batch (for conformity controlled cylinders).

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### FORTUNA® Graduated Cylinders

high form, **Class B**,  
borosilicate glass 3.3,  
hexagonal base

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**Packing unit:** 2 (\*1)

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### FORTUNA® Graduated Cylinders

high form, **Class A**,  
borosilicate glass 3.3,  
hexagonal base

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**Packing unit:** 2 (\*1)

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### FORTUNA® Graduated Cylinders


low form, **Class B**,  
borosilicate glass 3.3,  
hexagonal base

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### **FORTUNA®** Graduated Cylinders

high form, with PP-base,  
**Class B**, soda lime glass

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**Packing unit: 2**

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### FORTUNA® Mixing Cylinders

**Class A**, conformity approved,  
borosilicate glass 3.3,  
hexagonal base,  
with PP-stopper

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## **volac** Cylinders

All **volac** cylinders are made of borosilicate glass, have hexagonal feet and are printed in permanent amber stain or white ceramic glaze.

### Technical details:

- every spouted cylinder (graduated) is formed to ensure clean, non drip pouring
- every stoppered cylinder (mixing) has a tooled and ground neck fitted with a 'Polystop' PP interchangeable stopper
- all cylinders are individually calibrated and individually printed to ensure that the scale length exactly fits the calibration of each cylinder,
- the Class A cylinders are batch tested with a lot number, WAC cylinders are individually tested and numbered
- if a Works Certificate is required this will show the date of testing, test results at 5 points and the identification number of the cylinder.
- All Class A and B cylinders satisfy the requirements of International Standards (**EN ISO 4788, ASTM E1272-02**).

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
**volac**

### Graduated Cylinder

**Class B**, amber graduation,  
borosilicate glass,  
hexagonal base

**R250**

**volac**

### Graduated Cylinder

**Class B**, white graduation,  
borosilicate glass,  
hexagonal base

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### Graduated Cylinder

**Acc. to ASTM and ISO**

**Class A**, white graduation,  
borosilicate glass,  
hexagonal base

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
### **USP Graduated Cylinder**


**USP ASTM Class A,**  
white graduation,  
borosilicate glass,  
hexagonal base

**US263/WAC**  
with individual serial number  
and single certificate

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**volac**

### Mixing Cylinder

**Class B**, amber graduation, borosilicate glass, hexagonal base, with PP-stopper

R251

**volac**

### Mixing Cylinder

**Class B**, white graduation, borosilicate glass, hexagonal base, with PP-stopper

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**volac**

### Mixing Cylinder

**Acc. to ASTM and ISO  
Class A**, white graduation,  
borosilicate glass,  
hexagonal base,  
with PP-stopper

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### **volac** USP Mixing Cylinder

**USP ASTM Class A,**  
white graduation,  
borosilicate glass,  
hexagonal base,  
with PP-stopper

**US264/WAC**  
with individual serial number  
and single certificate

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**Packing unit:** 2 (\*1)

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### FORTUNA® Volumetric Flasks

FORTUNA® volumetric flasks 'Made in Germany' are produced to the highest quality standards in borosilicate glass only. They are manufactured in class A, conformity approved. Volumetric flasks are available with rim, with ground neck (TS), with PP- or glass stopper. (EN ISO 1042).

**Reading:** FORTUNA® volumetric flasks are calibrated to contain ('IN') and are offered with coloured ring mark. Reading of the volume is at the lowest point of the meniscus.

#### Technical Details:

- borosilicate glass 3.3
- accuracy class A
- with tolerance indications
- amber ring marks acid and alkali resistant
- indication of production batch

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### FORTUNA® Volumetric Flask

**Class A**, conformity approved,  
borosilicate glass 3.3

#### 1.500

amber graduation,  
with rim

#### 1.514-..-02F

blue graduation,  
with TS, without stopper

#### 1.514-..-04F

amber graduation,  
with TS, without stopper

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### **FORTUNA® Volumetric Flask**

**Class A**, conformity approved,  
borosilicate glass 3.3

#### **1.513-...-02F**

blue graduation,  
with glass stopper

#### **1.513-...-04F**

amber graduation,  
with glass stopper

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### **FORTUNA® Volumetric Flask**

**Class A**, amber stain diffused,  
conformity approved,  
borosilicate glass 3.3

#### **1.515**

white graduation,  
with PP-stopper

#### **1.516**

white graduation,  
with glass stopper

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### **FORTUNA® Volumetric Flask, Kohlrausch pattern**

**Class A**, without conformity approval, borosilicate glass

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### **volac** Volumetric Flasks

We offer a range of flasks from 1 ml to 5000 ml capacity, all manufactured to **Class A** (ISO and ASTM) standards, with permanent marking and tooled ground neck, fitted with interchangeable 'Polystop<sup>®</sup>' polypropylene stopper. All flasks are made of borosilicate glass and are available in clear or amber glass (**EN ISO 1042**).

#### **Class A Accuracy**

This range has been enlarged to cover all the former class B capacities. All flasks are individually tested for accuracy and are supplied with a batch test certificate and lot number.

#### **Class A Accuracy, Type WAC**

Each flask is individually tested and carries its own individual number, with works certificate of calibration.

#### **Class A according to USP**

Each flask is individually tested and carries its own individual number. A works certificate of calibration can be provided at extra cost.

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**volac**

### Volumetric Flask

**acc. to ISO and ASTM,  
Class A**, white graduation,  
individually calibrated, serially  
numbered, borosilicate glass,  
with Polystop®

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**volac**

### Volumetric Flask

**acc. to ISO and ASTM,  
Class A**, amber stain diffused,  
white graduation, individually  
calibrated, serially numbered,  
borosilicate glass, with  
Polystop®

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**RB258/AM**

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**RB258/AM/WAC**

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Capacity Inhalt	TS NS	RB258/AM	RB258/AM/WAC
		Art.-No. Art.-Nr.	Art.-No. Art.-Nr.
5 ml	10/19	RB258/AM/A/5	RB258/AM/WAC/A/5
10 ml	10/19	RB258/AM/B/5	RB258/AM/WAC/B/5
15 ml	10/19	RB258/AM/C/5	RB258/AM/WAC/C/5
20 ml	10/19	RB258/AM/D/5	RB258/AM/WAC/D/5
25 ml	10/19	RB258/AM/E/5	RB258/AM/WAC/E/5
50 ml	12/21	RB258/AM/F/5	RB258/AM/WAC/F/5
100 ml	12/21	RB258/AM/G/5	RB258/AM/WAC/G/5
200 ml	14/23	RB258/AM/H/2*	RB258/AM/WAC/H/2*
250 ml	14/23	RB258/AM/I/2*	RB258/AM/WAC/I/2*
500 ml	19/26	RB258/AM/J/2*	RB258/AM/WAC/J/2*
1000 ml	24/29	RB258/AM/K/2*	RB258/AM/WAC/K/2*
2000 ml	29/32	RB258/AM/L/5**	RB258/AM/WAC/L/5**

**Packing unit:** 5 (\*2, \*\*1)

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### **volac** USP Volumetric Flask

**USP ASTM Class A**, white graduation, each flask is individually tested and serially numbered and fully traceable satisfying USP/ASTM / ISO/BS/DIN accuracies, with Polystop®

#### **US258/WAC**

with individual serial number and single certificate

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### **volac** USP Volumetric Flask

**USP ASTM Class A**, amber stain diffused, white graduation, each flask is individually tested, serially numbered and fully traceable satisfying USP/ISO/BS/DIN accuracies, with Polystop®

**US258/AM/WAC**  
with individual serial number and single certificate

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## US258/AM/WAC

Capacity Inhalt	TS NS	Art.-No. Art.-Nr.
5 ml	10/19	US258/AM/WAC/A/5
10 ml	10/19	US258/AM/WAC/B/5
15 ml	10/19	US258/AM/WAC/C/5
20 ml	10/19	US258/AM/WAC/D/5
25 ml	10/19	US258/AM/WAC/E/5
50 ml	12/21	US258/AM/WAC/F/5
100 ml	12/21	US258/AM/WAC/G/5
200 ml	14/23	US258/AM/WAC/H/2*
250 ml	14/23	US258/AM/WAC/I/2*
500 ml	19/26	US258/AM/WAC/J/2*
1000 ml	24/29	US258/AM/WAC/K/2*
2000 ml	29/32	US258/AM/WAC/L/S**
5000 ml	34/35	US258/AM/WAC/M/S**

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volac

### Trapezoidal Flask

**Class A**, white graduation,  
borosilicate glass, with  
PP-stopper TS 7/16

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<sup>(1)</sup> not listed in EN ISO 1042 · nicht in EN ISO 1042

<b>Capacity</b> Inhalt	<b>TS</b> NS	<b>RB257</b>	<b>RB257/AM</b>
		<b>Art.-No.</b> Art.-Nr.	<b>Art.-No.</b> Art.-Nr.
1 ml	7/16	RB257/A/5	RB257/AM/A/5
2 ml	7/16	RB257/B/5	RB257/AM/B/5
3 ml	7/16	RB257/C/5 <sup>(1)</sup>	RB257/AM/C/5 <sup>(1)</sup>
4 ml	7/16	RB257/D/5 <sup>(1)</sup>	RB257/AM/D/5 <sup>(1)</sup>
5 ml	7/16	RB257/E/5	RB257/AM/E/5

<sup>(1)</sup> not listed in EN ISO 1042 · nicht in EN ISO 1042

**Packing unit: 5**

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### **FORTUNA® Precision Gas Syringes (Catheter Syringes)**

It is well known that gas volumetric analyses are indispensable to a better understanding of the fundamental principles of chemistry. Up to now, however, they have always been rather complicated as regards to the instruments involved (gas burette, Hempel's absorption pipette, explosion or slow burning gas pipette) and seen in the light of the difficulties facing the chemistry teacher, it was difficult to be sure of success with them.

The FORTUNA® precision gas syringe (catheter syringe) not only eliminates all these difficulties but makes tests so much quicker and easier to carry out so that it is already being used in both schools and universities in many countries. Its popularity is largely due to the simplicity of design and the exclusion of confining liquid, thus eliminating solution errors and enabling tests to be carried out with gases highly soluble in water.

#### **Technical Details:**

- clear Durobax® glass
- calibrated cylinders
- acid and alkali resistant amber stain diffused graduations
- FORTUNA® gas syringes are only for technical use

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2.012

2.000

The quantitative combustion of gases in the glow wire bulb is very much simplified too. Some of the tests which are important to teaching are as follows:

- analysis of inhaled and exhaled air (determination of breathing quotients)
- indirect water analysis
- combustion of elements for atomic weight analysis
- quantitative gas combustion ( $H_2$ ,  $CO$ ,  $NH_3$ , hydrocarbons)
- photosynthesis of hydrochloric acid

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### FORTUNA® Iodine number flask, Sendtner pattern

**2.120**

with hexagonal stopper,  
borosilicate glass 3.3

2.120

2.152

### FORTUNA® Graduated Cylinder for compressed volume

**2.152**

high shape, blue graduation,  
soda lime glass, round base

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**VTLB01- VTLB09**

### Density Bottle, acc. to Gay Lussac

pear-shaped, top of the stopper  
ground and polished,  
borosilicate glass 3.3 (DIN ISO 3507)

**2.500**  
uncalibrated

**2.510**  
calibrated to contain 1N (20°C),  
absolute capacity specified to  
 $\pm 0.001$  ml on each bottle

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
Gerber Milk Irish (1S66)	Gerber Milch Irish (1S66)	10.80 ml	R300/1
--------------------------	------------------------------	----------	--------

Gerber Acid	Gerber Säure	10.0 ml	R302
-------------	--------------	---------	------

Cream Pipette Double Bulb	Rahmpipette	5 g	R299
---------------------------	-------------	-----	------

Amyl Alcohol	Amylalkohol	1 ml	R304
--------------	-------------	------	------

Double Gerber pipette	Doppelte Gerberpipette	21.88 ml	R300/X2
-----------------------	---------------------------	----------	---------

Milk Pipette (BS755)	Milchpipette	17.6 ml	R300/17.6
----------------------	--------------	---------	-----------

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### **volac** Glass Pasteur Pipettes Made in the UK

All **volac** Glass Pasteur Pipettes are precision made from carefully selected tubing on fully automatic machines ensuring total uniformity of production. Tips and tapers are centred and can be checked by rolling the pipettes on a flat surface. Length and taper of the tips are carefully controlled. Constrictions for plugging are neatly formed at a uniform distance from the fire polished end. Constrictions are heat formed which avoids weaknesses or contamination of the glass found when moulding or tooling processes are used (ISO7712).

**volac** Glass Pasteur Pipettes are made with a slightly heavier wall thickness to provide the best quality product allied to greater strength to give lower breakages during transit and reduced risk of breakage in use or when fitting PVC teats. **volac** Glass Pasteur Pipettes are available either plain, plugged with non-absorbent roving or pre-plugged and sterilised by vacuum autoclaving at 134 °C.

#### Technical Details:

- available in standard sizes 150 mm, 230 mm, 270 mm, plain or plugged, sterilised or non-sterile
- '50 dropper' pipettes also available
- superior wall thickness as standard
- pasteurs 230 mm long can be supplied with heavy wall
- all standard pipettes can be supplied with 'sealed points' and break marks on request

#### Technical Details

##### D810, D812, D 812 HW (length 230 mm)

	Standard Past. Pip. mm	Heavy wall Past.Pip. mm
Top Body outer $\varnothing$	6.8 - 7.1	7.70 - 8.10
Wall thickness	0.5 - 0.53	0.80 - 0.90
Jet outer $\varnothing$	1.50 $\pm$ 0.2	

## Pasteurpipetten aus Glas Made in the UK

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#### D813

Teat for Pasteur Pipettes made of PVC, smooth for safety and durability

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#### D813/AD

Adaptor, for use with 1 ml disposable serological pipettes

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
**volac** Disposable Glass Serological Pipettes offer economic single use convenience leading to greatly reduced possibility of cross-contamination. Glass disposables have clear advantages over their plastic counterparts in areas of chemical and temperature resistance, accuracy, transparency and clarity of markings. All disposable pipettes are accurate to  $\pm 2\%$ .

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plugged, sterilised  
(in bulk packs)

#### D802/1/PS - D808/1/PS

plugged, sterilised  
(in sterile packs)

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### **Roving (available in three sizes wound on fibre spools)**

#### **R831**

for 1 ml pipettes:	approx. number of plugs produced 30,500, min. length 1150 m
for 2,5,10 ml pipettes:	approx. number of plugs produced 12,000, min. length 450 m
for wide top pip. and pasteur pipettes:	number of plugs produced 7,200, min. length 275 m

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Because we care.

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## 04. Syringes · Spritzen

## Classification of Medical Products

The development of European countries into a single market has made it necessary to follow common guidelines binding on all manufacturers. All customers are guaranteed that all medical products fulfil the quality and safety standards of the European Union. In the field of human medicine the Council of the European Union passed guidelines in 1993 called the MDD (Medical Directive Device). This law specifies that manufacturers have to evaluate and describe the risks and the proper use of medical products.

Poulten & Graf has established a quality management system according to DIN EN ISO 9001 et seq. and a further certification according to DIN EN ISO 13485.

These norms and specifications are updated and renewed continuously so please ask our medical products customer service for more information if required.

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## 04. Syringes · Spritzen

### FORTUNA® OPTIMA® All Glass Syringe

Medical device, CE,  
Durobax® glass, interchangeable

#### 7.102

clear glass, with glass tip,  
Luer, amber graduation

#### 7.140

clear glass, with Luer-Lock  
tip, amber graduation

#### 7.140-..-ATBG

amber stain diffused,  
with Luer-Lock tip,  
white graduation

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7.162

### FORTUNA® OPTIMA® All Glass Syringe, model ‚Tuberculin‘

Medical device, CE,  
Durobax® glass, long type,  
amber graduation

**7.162**

with glass tip, Luer

**7.168**

with Luer-Lock tip

### FORTUNA® OPTIMA® Glass and Metal Syringe

Durobax® glass, with steel  
piston, interchangeable,  
amber graduation, for  
technical use only

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

7.243

### Spare barrel for Glass and Metal Syringe FORTUNA® OPTIMA®

Durobax® glass, amber graduation,  
for technical use only

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## 04. Syringes · Spritzen

### FORTUNA® OPTIMA® Glass and Metal Syringe, model ‚Tuberculin‘

Durobax® glass, long type,  
with steel piston, interchangeable,  
amber graduation, for technical use  
only

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**7.264**

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### FORTUNA® OPTIMA® JANET, Washing Syringe

Durobax® glass, with metal piston,  
interchangeable, with 2 tips and  
metal holding bow, amber graduation,  
for technical use only

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## 04. Syringes · Spritzen

7.415

### HSW Single Use Syringes, 2-part, without needle, ISO 7886-1

sterile, CE

#### 7.410

with Luer tip

#### 7.415

with Luer-Lock tip

#### 7.418

with Catheter tip

### HSW Single Use Syringes, 2-part, with by-packed needle, ISO 7886 + 7864

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## 04. Syringes · Spritzen

7.452

### HSW Single Use Syringes, Insulin, without needle, ISO 8537

long model, with rubber gasket,  
sterile, CE

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**7.453**

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### Single Use Needles, ISO 7864

standard, sterile,  
non-corroding steel, CE

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**Packing unit:** 100**Verpackungseinheit:** 100

# Poulten & Graf Laboratory Products for your needs

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## 05. Laboratory Products · Laborprodukte

For all the following products we guarantee the quality of raw materials and precise performance of the finished products. The user of glass and plastic labware can be confident that they can count on laboratory equipment from Poulten & Graf.

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## Erlenmeyer Flask

with rim, white graduation,  
borosilicate glass

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### Erlenmeyer Flask with TS Neck, ISO 4797

**6.246/6.250**

white graduation,  
Borosilicate glass

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Borosilicate glass

**8.080**

with interchangeable plastic knob  
lid and porcelain plate

**8.100**

with tubes (outlet) and stopcock TS  
24/29 in lid and porcelain plate

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Borosilicate glass

**6.724**

cylindrical type, with conical,  
ground joint NS, PTFE-stopper,  
ungraduated

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## 05. Laboratory Products · Laborprodukte

### Glass stopper, DIN 12252

with hexagonal plate,  
hollow, borosilicate glass

**6.946**  
clear glass

**6.946-..B**  
amber glass

### PP-Stopper, DIN 12254

with sealing profile

**101444**

### Polystop®

interchangeable stopper,  
made of PE

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borosilicate glass 3.3

#### L180

tall form, with ground hollow  
blown and totally enclosed  
stopper

#### L181

tall form,  
with ground over cap

#### L182

low form, with ground hollow  
blown and totally enclosed  
stopper

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9.710

9.700

### Laboratory Bottle with screw thread, EN ISO 4796-1

#### 8.780

clear glass, white graduation, autoclavable, borosilicate glass, with cap

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### Amber Glass Bottle with screw thread

#### 9.700

non coated, 250 - 1000 ml: square 2500 ml: round

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#### 9.700

unbeschichtet, 250 - 1000 ml: vierkantig 2500 ml: rund

#### 9.710

plastic coated, 250 + 500 ml: square, 1000 + 2500 ml: round

#### 9.710

kunststoffbeschichtet, 250 + 500 ml: vierkantig, 1000 + 2500 ml: rund

	<b>8.780</b>
<b>Capacity</b>	<b>Art.-No.</b>
<b>Inhalt</b>	<b>Art.-Nr.</b>
100 ml	8.780-24
250 ml	8.780-36
500 ml	8.780-44
1000 ml	8.780-54
2000 ml	8.780-63

**Packing unit:** 10

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	<b>9.700</b>	<b>9.710</b>
<b>Capacity</b>	<b>Art.-No.</b>	<b>Art.-No.</b>
<b>Inhalt</b>	<b>Art.-Nr.</b>	<b>Art.-Nr.</b>
250 ml	9.700-02	9.710-02
500 ml	9.700-04	9.710-04
1000 ml	9.700-06	9.710-06
		9.710-07 <sup>(1)</sup>
2500 ml	9.700-08 <sup>(1)</sup>	9.710-08 <sup>(1)</sup>

<sup>(1)</sup>round

<sup>(1)</sup>rund

**Packing unit:** 1

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## 05. Laboratory Products · Laborprodukte

RP705

RP705/AM

GS705

GS705/AM

### Plain Narrow Mouth Reagent Bottle

#### RP705

clear glass, with Polystop®

#### RP705/AM

amber glass, with Polystop®

#### GS705

clear glass, with solid amber  
glass stopper

#### GS705/AM

amber glass, with solid amber  
glass stopper

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Capacity Inhalt	TS NS	RP705	
		Art.-No. Art.-Nr.	Art.-No. Art.-No.
30 ml	16/16	RP705/A/10	-
50 ml	16/16	RP705/B/10	RP705/AM/B/10
100 ml	16/16	RP705/C/10	RP705/AM/C/10
250 ml	24/20	RP705/D/10	RP705/AM/D/10
500 ml	24/20	RP705/E/2*	RP705/AM/E/2*
1000 ml	29/32	RP705/F/2*	RP705/AM/F/2*

Capacity Inhalt	TS NS	GS705	
		Art.-No. Art.-Nr.	Art.-No. Art.-No.
30 ml	16/16	GS705/A/10	-
50 ml	16/16	GS705/B/10	GS705/AM/B/10
100 ml	16/16	GS705/C/10	GS705/AM/C/10
250 ml	24/20	GS705/D/10	GS705/AM/D/10
500 ml	24/20	GS705/E/2*	GS705/AM/E/2*
1000 ml	29/32	GS705/F/2*	GS705/AM/F/2*

**Packing unit:** 10 (\*2)

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## 05. Laboratory Products · Laborprodukte

SCG705

### Safety Coated Surfisil® Reagent Bottle

made of soda lime glass

#### SCP705

clear glass, with Polystop®

#### SCP705/AM

amber glass, with Polystop®

#### SCG705

clear glass, with solid amber glass stopper

#### SCG705/AM

amber glass, with solid amber glass stopper

### Indicator Pipette Dropping Bottle, 50 ml

#### SL920

clear glass, available labelled with 6 standard titles

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<b>Capacity</b> Inhalt	<b>TS</b> NS	<b>SCG705</b>		<b>SCG705/AM</b>	
		<b>Art.-No.</b> Art.-Nr.		<b>Art.-No.</b> Art.-No.	
30 ml	16/16	SCG705/A/10		-	
50 ml	16/16	SCG705/B/10		SCG705/AM/B/10	
100 ml	16/16	SCG705/C/10		SCG705/AM/C/10	
250 ml	24/20	SCG705/D/10		SCG705/AM/D/10	
500 ml	24/20	SCG705/E/2*		SCG705/AM/E/2*	
1000 ml	29/32	SCG705/F/2*		SCG705/AM/F/2*	

<b>Title</b> Aufschrift	<b>SL920</b>	
	<b>pH-Range</b> pH-Bereich	<b>Art.-No.</b> Art.-No.
Litmus	5 - 8	SL920/B/L28
Methyl Orange	2.9 - 4.6	SL920/B/L29
Methyl Orange Screened	2.9 - 4.6	SL920/B/L30
Phenolphthalein	8.3 - 10	SL920/B/L31
Starch Indicator	-	SL920/B/L32
Universal Indicator	3 - 11	SL920/B/L33

## 05. Laboratory Products · Laborprodukte

### Laboratory Bottle, labelled, safety coated

#### SL921

clear glass with Polystop®

- A indicates amber glass bottle
- G indicates amber glass bottle with amber glass stopper

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<b>Title</b>	<b>Art.-No. 250 ml</b>	<b>Art.-No. 500 ml</b>
<b>Aufschrift</b>	<b>Art.-Nr. 250 ml</b>	<b>Art.-Nr. 500 ml</b>
Acetic Acid Dil	SL921/D/L01	SL921/E/L01
Acetic Acid, Glacial	SL921/D/L02	SL921/E/L02
Ammonia Soln. Dil	SL921/D/L03	SL921/E/L03
Ammonia Conc.	SL921/D/L04	SL921/E/L04
Barium Chloride	SL921/D/L05	SL921/E/L05
Bromine Water (A,G)	SL921/D/L06	SL921/E/L06
Copper (II) Sulphate	SL921/D/L07	SL921/E/L07
Ethanol	SL921/D/L08	SL921/E/L08
Fehlings Soln. No. 1	SL921/D/L09	SL921/E/L09
Fehlings Soln. No. 2	SL921/D/L10	SL921/E/L10
Hydrochloric Acid Dil.	SL921/D/L11	SL921/E/L11
Hydrochloric Acid Conc. (G)	SL921/D/L12	SL921/E/L12
Lead (II) Acetate	SL921/D/L13	SL921/E/L13
Lime Water	SL921/D/L14	SL921/E/L14
Nitric Acid Dil.	SL921/D/L15	SL921/E/L15
Nitric Acid Conc. (A,G)	SL921/D/L16	SL921/E/L16
Potassium Dichromate (VI)	SL921/D/L17	SL921/E/L17
Potassium Hydroxide	SL921/D/L18	SL921/E/L18
Potassium Iodide	SL921/D/L19	SL921/E/L19
Potassium Manganate (VII)	SL921/D/L20	SL921/E/L20
Propanone (Acetone)	SL921/D/L21	SL921/E/L21
Silver Nitrate (A)	SL921/D/L22	SL921/E/L22
Sodium Carbonate	SL921/D/L23	SL921/E/L23
Sodium Hydroxide Dil.	SL921/D/L24	SL921/E/L24
Sodium Hydroxide Conc.	SL921/D/L25	SL921/E/L25
Sulphuric Acid, Dil.	SL921/D/L26	SL921/E/L26
Sulphuric Acid Conc. (G)	SL921/D/L27	SL921/E/L27

## 05. Laboratory Products · Laborprodukte

L188/AM

### Polystop® Dropping Bottle

with Polystop®, with 3 piece stopper (including a strong ungraduated glass pipette, a durable vinyl teat and a PP-stopper)

**L188**

clear glass

**L188/AM**

amber glass

### Dispensing Bottle

with 3 piece polystopper (including a strong amber graduated dropping pipette, a durable vinyl bulb and a PP-stopper)

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### Balsam Bottle

#### L669

clear glass, with loose glass cap and glass rod

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### Cedar Wood Oil Bottle

#### L670

clear glass, has a C 16 Polystopper fitted with a glass rod

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Biological Oxygen Demand

#### GS748

clear glass, has glass peg stopper adjusted to leave a 'well' for water seal and a pointed cone to prevent air entrapment

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**GS748/ST**

for 250/500 ml

#### **Polystop® 3 piece stopper for dropping bottles**

**L188/SAV**

including a strong ungraduated glass pipette, a durable vinyl teat and a PP-stopper

#### **Plain Pipette**

**L188/GP**

for before mentioned assembly

#### **3 piece polystopper for dispensing bottles**

**L700/SAV**

including a strong amber graduated dropping pipette, a durable vinyl bulb and a PP-stopper

#### **Plain Pipette for Dispensing Bottle**

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
### Flat hat stopper, amber glass GS705/././ST

### Storage Bottle, PTFE 10.1305

with screw cap (suitable for  
FORTUNA UNIVERSAL® HF dispenser)

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## 05. Laboratory Products · Laborprodukte

10.1290-..B

### Spare Bottle for Burettes Dr. Schilling pattern

polyethylene, without screw cap

**10.1290**

white

**10.1290-..B**

amber

### Plastic Transfer Pipettes

Plastic Transfer Pipettes are precision products made from high quality, chemically inert, non-toxic low-density polyethylene, ideal for handling liquids safely in controlled, reproducible drop sizes in all types of laboratory tests and analyses.

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**Poulten & Graf Ltd.**

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**What is understood by quality?**

In general usage the term 'Quality' has a number of meanings. Glass products of first class quality might be called this because the finish is good, the raw material used is good, the packaging is good, the appearance is good or the accuracy is good. Quality is defined as having all those characteristics and features of a product, which refers to its suitability for the fulfilment of given requirements, that is 'fit for purpose'. So Quality is related to the need and the specification set.

**The Quality Management of Poulten & Graf**

Product quality is directly dependent on the quality of process. Quality assurance must begin in the very first phase of the product creation process and cover all following phases. This global view of quality assurance, involving all aspects of production, design, testing, packing and conformity with both internal standards and externally set specifications, is the basis of our quality management. This is reflected in the approvals of DIN EN ISO 9001/2008, which we have achieved and which we maintain.

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# 06. Technical Support · Technische Details

## Certification or attestation of conformance

Since January 01, 2015, new German Metrology and Verification Regulations have come into force. They replace the previous German Weights and Measures Regulations. For Class A and Class AS Volumetric Glassware which are produced in Germany, a new conformity sign is being used. The H-sign is now replaced by the **DE-M**-sign.

As before, the manufacturer is allowed to certify the conformity on basis of the PTB-regulations (National Metrology Institute). Government institutions will ensure that the regulations are being followed. Therefore, nothing has changed in the production and testing methods. If, required, a manufacturer's declaration of conformity can be enclosed with class A and AS volumetric products.

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## Works Certification of Accuracy of Volumetric Glassware

All  Class A glassware is always individually tested. Test records are kept in house with the manufacturing records. However individual Works Certificates of Accuracy can be provided for each piece if required; each piece of glass has its own individual number engraved on it and the Certificate supplied will have the corresponding number, the date and the relevant test results.

Such Works Certificates of Accuracy (look for the designation WAC in the product codes) should not be confused with Certificates of Conformity which refer to batch testing only.

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## The system for lot numbers on FORTUNA® volumetric glassware

All FORTUNA® conformity approved volumetric glassware is delivered with the FORTUNA® Batch Certificate of Performance. A regular control of all testing equipment is required by DIN EN ISO 9001. The requirement for traceability is met by the batch certificate and with the easily identifiable lot numbers (e.g. LOT B008) on FORTUNA® volumetric glassware, which is produced in Wertheim.

Beside the conformity sign all FORTUNA® volumetric glassware is marked with an individual lot number that allows an exact classification back to the batch certificate.

If your certificate is lost, you may ask for a copy by informing us of the product details and the batch number or request it online.

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R9009

To define the different classes (Class A, Class AS, Class B, etc.) for glass or 'performance' for pipetters and dispensers the following terms are used: 'accuracy' and 'coefficient of variation' (precision). Sometimes we also speak about 'tolerances' which are specified in the ISO-standards. This means the permitted deviation of the measured values from the nominal volume. By accuracy we mean corresponding to the average value of a series of measurements adjusted for temperature and pressure with the desired value. Precision is a measure of the variation from the nominal value achieved through a series of measurements at one setting. The coefficient of variation is the % value of the precision.

Measurement of volumes is part of the routine work in a laboratory. The most common volumetric glassware such as bulb- and graduated pipettes, burettes, volumetric flasks and graduated cylinders form part of the basic equipment in a laboratory. Inscriptions and graduations on volumetric glassware should be easy to read and acid- and alkali-resistant. The complete range of FORTUNA and  volumetric glassware fulfils these requirements.

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
### Calibration of Volumetric Glassware

- **EX** calibrated to deliver. The delivered liquid corresponds with the volume indication. The liquid that remains on the wall or in the tip of the glassware is taken into consideration.
- **EX + s** calibrated to deliver after a waiting time (for pipettes 5 seconds, for burettes 30 seconds). Keeping to the waiting time is important to avoid mistakes in measurement.
- **IN** calibrated to contain. The contained liquid corresponds with the volume indication. Example: flasks and cylinders.
- **20°C** The calibration has been done at 20°C. Volumetric glassware changes volume about 0,1 % each degree Celsius.

### Precision Classes

- **A** tolerance according to ISO standard. Class A is the most precise class. Glassware of class A is suitable for official calibration, this means it can be tested by either the German or the UK Weights and Measures Offices. This is useful for control of measuring instruments according to DIN ISO 9001.
- **AS** as class A, but for pipettes and burettes with fast delivery. The waiting time is much shorter than for class A.
- **B** the tolerances of class B are about twice those of class A and AS. Volumetric glassware of class B is suitable for routine work in the laboratory.

### Meniscus

Value is read at the lowest point of the liquid surface, which means when the lowest point of the meniscus touches the upper line of the graduation mark (see picture 1).

For burettes with Schellbach stripes the reading has to be made at the level where the two wedge-shaped points are touching (see picture 2). When reading the scales it is important that one's eyes are at height of the fluid level. Otherwise it results in reading errors (parallax errors).

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
### A Guide to calibrate Micro Pipettors and Bottle Top Dispensers (compiled from BS ISO 8655: Part 1-6)

#### Test Lab Requirement

- Lab should ideally be air conditioned.
- Temperature and air pressure should be stable.
- Air movement should be minimised.
- Weighing/balance surface should be free of vibration.

#### Test Lab Equipment

- Air pressure barometer
- Certified Class A thermometer
- Accurate balance (5 decimal places)  
properly calibrated to National Standards
- Distilled water at stable room temperature

#### Pipetting Practice for Accurate Results

Use high quality tips preferably of a brand to match the pipettor. Ensure tips do not leak (hold pipettor vertically, tip filled with liquid, without dripping for at least 20 secs). Rinse and pre-wet each tip with distilled water at least 3 times before use. Fill and dispense slowly and smoothly. Fill with pipettor tip not more than 3-4 mm below liquid surface and pipettor held vertical (fig. 1) Dispense with pipettor tip against wall of container and pipettor held at 45° to the vertical and draw tip up the container wall one or two seconds after second blow out (fig. 2).

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

## 06. Technical Support · Technische Details

### Pipetting - Balance Readings

Tare a container on the balance. Using a wetted tip, pipette/dispense into the tared container on the balance. Read weight of dispensed water. Read temperature of distilled water. Read the air pressure. For setting calibration test 4 times at each volume, for checking/certifying calibration test minimum 10 times and for testing coefficient of variation test 30 times.

### Calculation of Accuracy

- Calculate mean weight  $\bar{W}$  from the number  $n$  of individual weighings  $W_i$  summed to give  $\sum W_i$  using the equation:

$$\bar{W} = \frac{\sum W_i}{n}$$

- Calculate mean volume  $\bar{V}_t$  of liquid samples at test temperature  $t$ :

$$\bar{V}_t = \bar{W} \cdot Z$$

where  $Z$  is the conversion factor for the density of water at given temperatures and air pressure - see table (page 294)

N.B. If an electronic balance is used, the effects of evaporation are minimal and can be ignored; if a mechanical balance is used mean evaporation loss  $\bar{e}$  needs to be estimated and added to  $\bar{W}$ .

- Calculate inaccuracy in percent of the mean  $\bar{A}_t$  at test temperature  $t$  using  $\bar{V}_t$  mean volume and  $V_o$  nominal volume by equation:

$$A_t = \bar{V}_t - V_o \cdot \frac{100}{V_o}$$

### Index of Symbols

$\bar{W}$  mean weight

$W_i$  results of individual weighings

$n$  number of weighings

$t$  water temp. of test

$\bar{V}_t$  mean volume at test temperature

$V_o$  nominal volume

$\bar{e}$  mean evaporation loss

$Z$  Conversion factor (see table)

$\bar{A}_t$  Accuracy (inacc.) as percent of mean

$S$  Standard deviation of weighings (or volumes)

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# 06. Technical Support · Technische Details

## Calculation of Precision (Repeatability) - Coefficient of Variation

- Calculate standard deviation **S** by the equation:

$$S = \sqrt{\frac{\sum (W_i - \bar{W})^2}{n-1}} \quad (\text{n should be not less than 10 and should ideally be 30})$$

- Calculate coefficient of variation **CV**:

$$CV = \frac{S}{\bar{W}} \cdot 100$$

If the balance is attached to computing equipment  $\bar{W}_i$  and **W** may be replaced by equivalent volumes.

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## Calibration of Volumetric Glassware (see BS 6696)

### Test Lab Requirement

- Full temperature control (air conditioning)
- Atmospheric Pressure Barometer
- Certified Thermometer
- Accurate calibrated weighing balance

### Cleaning and Drying

Wash, clean and dry glassware thoroughly. Ensure that a satisfactory method of cleaning is employed, since any contamination will effect the meniscus shape/level.

During calibration it is critical that all internal surfaces are clean and dry to avoid moisture droplets or contaminants.

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
### Timing Delivery

Delivery glassware should be timed before calibration, to ensure that total delivery time is within tolerances specified by the relevant BS and ISO standards.

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### Weigh and Tare

With delivery glassware such as pipettes and burettes, weigh and tare a suitable sized weighing vessel.


Glassware should be maintained at a controlled room temperature (air conditioning).

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
Fill with purified water a few millimeters above the graduation mark. Pipettes and burettes should be filled from below, i.e. through the tip. If burettes are filled from the top for added convenience, take care to avoid wetting the walls above the upper graduation mark and check that no air is trapped in the stopcock. For pipettes and burettes control the flow of excess liquid from the tip to allow the meniscus to settle exactly on the graduation mark. For flasks and cylinders remove excess water to match meniscus and graduation mark exactly using a pipette and bulb.

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
When setting the meniscus ensure that the glassware is vertical. Ensure parallax error is avoided by viewing the meniscus from an eyelevel height against a white background.

In addition a black shade should be secured either around the article to be tested or placed behind it. The shade should be below the setting height (approx 1 mm lower). Ensure the meniscus edge is smooth, straight and level.

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## 06. Technical Support · Technische Details

For graduated/bulb pipettes and burettes, dispense water into a tared weighing vessel and weigh. For flasks and cylinders weigh the tared and filled glassware. Record the result from the balance. Use a certified Class A thermometer to measure the water temperature. Add the water temperature correction figure to the net weight obtained.

Measure the barometer reading and air temperature and add (or subtract if minus) the correction figure. Please ask us for the tables.

This produces the final gravimetric result which has already been converted to the volumetric figure in millilitres.

All of the procedure described must be followed to ensure that consistently accurate results are obtained.

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## Environment and safety

Current discussions on environmental protection, dangerous materials etc. have made laboratory technicians more sensitive in the selection of laboratory instruments. Safety aspects remain the primary concern while handling dangerous or hazardous materials. Laboratory equipment must be extremely safe and reliable. All regulations must be observed to ensure protection of both humans and the environment.

We offer a repair service for a number of products:

- bottle top dispenser POLYFIX®, OPTIFIX®,
- electronic burette OPTIFIX® titrier,
- automatic dosing system OPTIMAT® 2,
- pipettors.

Please remember that when returning any instrument for repair it must be accompanied by a written certificate of decontamination.

When disposal of waste is involved we ask you to consider the necessary regulations. These apply also to disposal of single use products.

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## 06. Technical Support · Technische Details

### Important safety instructions:

1. Laboratory instruments and equipment should be checked on purchase as being in good condition and working order before use.
2. Laboratory instruments should be checked for damage as part of a daily routine to avoid risk of injuries such as cuts, burns and infections.
3. Laboratory glassware which shows any sign of damage should be disposed of in a safe manner.
4. Sudden changes in temperature should be avoided while handling glassware.
5. Glass equipment should be inspected thoroughly before using in a vacuum or under pressure.
6. Sudden changes in pressure should be avoided. Vessels should be vented smoothly when used in a vacuum. Flat bottomed glassware should not be used under pressure or in a vacuum.

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## 06. Technical Support · Technische Details

### Cleaning of Laboratory Glassware

Laboratory glassware should be washed before first use.

#### Cleaning by hand

The common cleaning method is to wipe and rub the glassware with a cloth or sponge soaked in a cleaning solution. Abrasive cleaners and sponges should not be used because they can damage the surface of the glass.

#### Machine cleaning

Cleaning in a machine is more gentle for glassware than soaking it. Glass comes into contact with the dishwashing liquid for relatively short periods when it is being sprayed onto the surface of the glass.

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## Cleaning of Bottle Top Dispenser

To ensure proper operation and life of the OPTIFIX® dispensers, we recommend the following cleaning procedure to provide you with years of trouble free service.

- Prepare a standard cleaning solution and hot water.
- Remove the piston from the Cylinder via the volume setting knob. Remove discharge tube, intake tube and adapter.
- Put all parts into the cleaning solution for 2-4 hours or overnight to allow a good soaking.
- After soaking, use a clean glassware brush to remove any further dirt if necessary.
- Rinse all dispenser parts with distilled water and air dry.
- Reassemble the unit and put the OPTIFIX dispenser back into service.

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---

### Disinfection

Laboratory equipment that has been in contact with infectious material should be disinfected before use, for the protection of laboratory workers. It can be washed by hand with a cleaner-disinfectant. Further it can be cleaned by using physical or thermal methods (for 10 minutes at 93°C as required by the BGA). If necessary the glassware can be steam sterilized.

### Steam sterilization

Steam sterilization is 'the destruction or irreversible inactivation of all microorganisms of 120°C and 2 bar'. At a sterilization temperature of 121°C the glassware should be sterilized for 20 minutes minimum effective application time. Laboratory equipment should be always cleaned carefully before steam sterilization, otherwise any soiling will bake on during the steam sterilization. If there are any chemicals on the surface of the glassware during the sterilization they can damage the glass due to the high temperatures. Containers should be opened during the sterilization to avoid a build-up of pressure. The steam has to have unrestricted access to all contaminated points to ensure effective steam sterilization.

### Raw materials we use

**Soda Lime glass** finds its use in the foodstuffs industry as bottle and preserve glass as well as in the laboratory industry. The chemical composition of this glass consists of 71-75 % sand ( $\text{SiO}_2$ ), 12-16 % soda and 10-15 % lime. Soda lime glass is subject to acute stress under heat shock and has a relatively high coefficient of expansion. It is particularly resistant to alkali liquids.

**Borosilicate glass** differs from soda lime glass by using a higher proportion of sand ( $\text{SiO}_2$ ). The composition of this glass contains 70-80% silicon dioxide ( $\text{SiO}_2$ ), 7-13% boron trioxide ( $\text{B}_2\text{O}_3$ ), 4-8%  $\text{Na}_2\text{O}$  and  $\text{K}_2\text{O}$ , as well as 2-7% alumina oxide ( $\text{Al}_2\text{O}_3$ ). Glass with such composition is characterised by a high resistance against acidic attack, low thermal expansion and excellent temperature stability. This glass is particularly suitable for laboratories and chemical use where contact with aggressive materials is likely to be prolonged.

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## 06. Technical Support · Technische Details

The POLYFIX® Dispenser is made of Polypropylene and glass.

The OPTIFIX® dispenser is made of PTFE and glass, while the SMART pipettor is made of Polypropylene and Polycarbonate.

### PP - Polypropylene

When polymerizing, propylene develops a product which contains methyl side groups of monomers of the paraffin chain. The volume requirement due to the regular methyl side chains is relatively large. For this reason polypropylene is the plastic with the lowest specific weight in compact form. It is similar to polyethylene, and may be used instead of polyethylene, where it particularly depends on good heat resistance, high impact strength, inherent stability during stress or particularly low specific weight. Polypropylene also has a good resistance to many inorganic reagents as well as the less aggressive organic solvents.

### PTFE - Polytetrafluorethylene

PTFE has no thermoplastic characteristics although it consists of linear, chemical chain molecules which are not linked. It is not possible to handle or form PTFE in the usual thermoplastic way. Nevertheless PTFE is a thermoplast because it can be sintered. For the production of shaped parts from PTFE the polymer is cold pressed as powder under high pressure into the desired form and sintered at temperatures around 380 °C. PTFE is attacked only by elemental fluorine and chlorine tri fluoride at higher temperatures and pressures as well as by melting alkali metals. Otherwise it is resistant to all chemicals. A further advantage of polytetrafluorethylene: it has the lowest coefficient of friction of all solids. However above 350 °C PTFE burns in air to form noxious gases.

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# Poulten & Graf Index, Terms and Conditions

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# 07. Index/Terms & Conditions · Register/AGBs

<b>A</b>			
Accessories for OPTIMAT®	56-59	Dispensing pumps for OPTIMAT®	56+57
Adaptors, taper size, for OPTIFIX®	46+47	Disposable glass pipettes, VOLAC	224+225
Adaptors, thread GL, for OPTIFIX®	46+47	Drain tube for dispensers	48+49
Adaptors for Pasteur Pipettes	222+223	Dropping bottles, glass	268-271
Addresses	5	Dr. Schilling burettes	102+103
All glass syringes	230-233	<b>E</b>	
Amber glass syringes	230+231	ECO microliter pipettes, VOLAC	70+71
ASTM glassware 170-175,180+181,196-201		Electronic bottle top titrator titrier	49-55
Automatic burettes	98-103	Erlenmeyer flasks, glass	246-249
<b>B</b>		<b>F</b>	
Balsam bottles	270+271	Filter Flasks	252+253
Bang burettes	96+97	Flasks, Erlenmeyer,glass,grad.	246-249
BASIC-Dispenser	30+31	Flasks, Kohlrausch	194+195
Beakers, FORTUNA®, glass	250+251	Flasks, volumetric, FORTUNA®	184-195
Blow-out pipettes 134+135+216+217		Flasks, volumetric, VOLAC	196-207
BOD bottles	270+271	Foot manipulator, f. OPTIMAT®	58+59
Bottles for disp., made of PTFE	274+275	Funnels, glass	248+249
Bottles made of PP	276+277	<b>G</b>	
Bottles with thread	260+261	Gas syringes, FORTUNA®	208-211
Bottle top dispenser OPTIFIX®	26-45	Glass stoppers 256+257+272+273	
Bottle top dispenser POLYFIX®	21-25	Glass syringes	230-233
Bulb pipettes, FORTUNA®	144-151	Glass and metal syringes 232+233+236+237	
Bulb pipettes, VOLAC	152-157	Glass Material	304+305
Burette stopcocks	122+123	Grad. cylinders, FORTUNA®	158-169
Burettes, FORTUNA®	90-103	Grad. cylinders, VOLAC	170-183
Burettes, VOLAC	104-123	Grad. pipettes, FORTUNA®	124-131
Butyrometer, VOLAC	214-215	Grad. pipettes, VOLAC	132-141
<b>C</b>		Grad. pipettes, VOLAC, disposable	224+225
Calibration guide + Service	288-299	<b>H</b>	
Carousel for Microliter Pipettes	68-69	Humidity absorber for OPTIFIX®	46+47
Catheter syringes	208-211	<b>I</b>	
Cleaning instructions	302+303	Injection needles	240+241
Cedar wood oil bottles	270+271	Intake tubes for OPTIFIX®/POLYFIX®	42+43
Conditions of Sale	316-329	Intake tubes for pumps OPTIMAT® 2	58+59
Cylinders for compressed volume	212+213	Iodine number flasks Sendtner pat.	212+213
Cylinders, graduated, FORTUNA®	158-169	Indicator pipette dropping bottles	264+265
Cylinders, graduated, VOLAC	170-183	<b>J</b>	
<b>D</b>		JANET wash syringes	236+237
Daffert burettes	120+121	<b>K</b>	
Density bottles	214+215	Kohlrausch flasks	194+195
Desalting unit LAB-IoN	76+77	<b>L</b>	
Desiccators	252+253	LAB-ION L 2 desalting unit	76+77
Discharge tubes for OPTIFIX®	44+45	Labelled bottles, glass	266+267
Discharge tube sets for OPTIMAT® 2	58+59	Laboratory bottles, ISO-thread, glass	260+261
Dispensers FORTUNA®OPTIFIX®	26-45		
Dispensers FORTUNA® POLYFIX®	21-25		
Dispensing bottles	268+269		
Dispensing burettes	120+121		

Laboratory bottles	260-271		
		M	
Micropipettors VOLAC SMART, ECO	60-71		
Milk pipettes	216+217		
MiniPipet, fixed volume micropipettor	70+71		
Mission Statements	8+9		
Mixing cylinders, FORTUNA®	168+169		
Mixing cylinders, VOLAC	178-183		
Multichannel pipette, VOLAC	66-69		
		N	
Needles	240+241		
Nessler tubes, VOLAC	216+217		
		O	
OPTIFIX® dispenser	26-45		
OPTIMA® syringes	228-237		
OPTIFIX® titrier	49-55		
Original FORTUNA® graduated pip.	130-132		
Original FORTUNA® bulb pipettes	150+151		
Overflow pipettes Daffert pattern	120+121		
		P	
Pasteur pipettes, glass, VOLAC	218-223		
Pellet automatic burettes	98-101		
Petri Dishes	254+255		
Pipettes	124-157		
Pipettes, single use, glass	224+225+274+275		
Pipette tips	72+73		
Pipetting aids	74+75+142+143		
Pipette dropping bottles	260-271		
Pipette fillers	74+75+142+143		
PI-PUMP	142+143		
Plastic Transferpipettes	276+277		
POLYFIX® dispenser	21-25		
Polystoppers	256+257		
PTFE-bottles for dispensers	274+275		
Pumps for OPTIMAT®	56+57		
PVC teats	222+223+274+275		
PVDF/PTFE-replacement stopcocks	122+123		
		Q	
Quality	10+11+279-287		
		R	
Reagent bottles VOLAC	262-265		
Record syringes	232-237		
Reductase tubes	216+217		
Replace cartridge for LAB-ION L2	76+77		
Reservoir bottles f. burettes	100+101+276+277		
Roving for pipettes	226		
		S	
SAFETY/SAFETY S-dipsenser	36-39		
Safety instructions	300+301		
Safety stopcock	44+45		
Serological pipettes, glass	224+225		
Serial cable for OPTIMAT®	58+59		
Single use needles	240+241		
Single use syringes	238+239		
SMART Micropipettor	60-69		
SOLVENT dispenser	32+33		
Spare barr. for glass/metal syringes	234+235		
Spare barrels for JANET syringes	236+237		
Spare cylinders for OPTIFIX®	40+41		
Spare parts for dropping bottles	272+273		
Spare parts for OPTIFIX®	40-45		
Spare parts for POLYFIX®	24+25		
Spare parts for Titrier	54+55		
Spare parts for water distillator	78+79		
Spare reservoir f. automatic burettes	100+101		
Stopper caps for OPTIFIX®	44+45		
Stoppers	256+257+274+275		
Syringes	228-241		
		T	
Tissue culture pipettes	136+177		
Teats	222+223		
Tips for Micropipettes	72+73		
Titratör OPTIFIX® titrier	49-55		
Trademarks	12		
Trapezoidal flasks	206+207		
Tuberculin glass syringes	232+233+236+237		
Tubing, PTFE	42+43		
		U	
Univ.drain tube for dispensers	48+49		
Universal HF dispenser	34+35		
USP volumetric glassware	112+113+140+141+156 +157+176+177+182+183+202+203+204+205		
		V	
Volume setting system for OPTIFIX®	42+43		
Volumetric flasks, FORTUNA®	184-195		
Volumetric flasks, VOLAC	196-207		
Volumetric glassware	81-226		
		W	
Wash Syringes, JANET	236+237		
Water distillator	78+79		
Weighing bottles, glass, VOLAC	258+259		
Wanderer, universal drain tube	48+49		

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# 07. Index/Terms & Conditions · Register/AGBs

## **Poulten & Graf Limited – Conditions of Sale**

### 1. General conditions

Our General Terms of Business are binding for all deliveries and services. We will not accept general terms of business of our customers, which differ from our terms of sale. Additional terms have to be confirmed by us in writing. Only English Law and English Text will be applicable. Place of performance is Barking. The place for jurisdiction, also for bill transactions, is England. All important data necessary for business transactions may be stored on computer. All information, which is not public knowledge and which results from the business relationship with our company, must be kept secret from third parties.

### 2. Offers

Our offers are always for information. In order to be valid, all transactions, orders and other agreements require our confirmation in writing. Our written order confirmation or agreement shall determine the nature and scope of the delivery. We reserve the right to make technical alterations to the offered products.

### 3. Risk, conditions of delivery and non-acceptance of delivery

The risk of ownership will be transferred to the buyer, as soon as the goods leave our company even in the event of delivery carriage-paid. Complaints can only be taken into account if submitted within 8 days of receipt of the goods. Delivered goods will only be accepted for return in their original packing. Illustrations and brochures etc. are not binding. Measurements and utilisation data are only to be regarded as approximate and are not binding.

Deliveries will be made ex factory, packing and insurance excluded, if no other agreement exists. Packing material of any kind may be returned to us within the range of legal regulations applicable to us. We will not pay for the cost of return of packing materials. In the event that dispatch or delivery is delayed at the request of the customer or for reasons for which he is responsible, the risk shall pass to the customer for the period of the delay, however at the request of the customer, we will arrange insurance cover that he requires at his expense.

Partial deliveries are permissible. Surplus quantities within the normal range (up to 15 %) are permissible and must be accepted and paid for by the customer. In case of delayed acceptance of the goods, we may sell the items at the customer's expense to a third party or we may ask for compensation.

### 4. Delivery times and Unforeseen Circumstances

Indicated delivery times are for information only, unless a binding delivery time has been fixed in writing. Deliveries may be delayed due to Act of God or strikes, interruption of work, shortage of raw materials or failure to deliver by our suppliers. Liability for delayed delivery will be limited to proven, gross negligence, caused by us.

### 5. Special productions

Products, which are manufactured according to drawings, samples or instructions of the customer, cannot be returned. This includes in particular individual software for computers,

delivered by us.

All kinds of software delivered are subject to copyright and may be neither copied, modified nor presented in public without our prior written permission. Additional quantities within the normal range (up to 15 %) are permissible and must be accepted and paid for by the customer, unless exact quantities have been agreed upon in writing. The customer is responsible for ensuring that products manufactured according to his instructions, do not interfere with any trade right, copyright or other rights of a third party. The customer must compensate the company for all damages resulting from any violation of such rights.

## 6. Payment conditions, prices, delayed payment

Products will be invoiced at the prices in effect on the day of delivery. Goods for delivery in UK will be in pounds sterling. All other prices will in Euros Prices are ex factory excluding packaging and transportation insurance with immediate payment, unless otherwise provided.

We reserve the right to include additional payment conditions for new customers. The Price for the Products must be paid in full and without setoff In case of delayed payment, we reserve the right to charge the customer with all costs as well as interest at 3 % above the Allied Irish Bank base rate from time to time.

## 7. Loss or Damage in Transit

We shall only accept responsibility for:

7.1. Damage to the goods caused in transit if the same is externally visible and is notified to us and the carrier (if not delivered by us) within 3 days of receipt of the goods by the customer.

7.2. An actual or apparent discrepancy between any delivery note and items delivered if the same is notified to us and the carrier (if not delivered by us) within 3 days of receipt of the goods by the customer.

7.3. Non arrival, if the same is notified to us within 3 days of any specific delivery date notified to the customer by us. If no specific delivery date is notified by us to the customer, then the customer must give notice within 21 days of the date of any invoice for the goods sent by us to the customer.

7.4. Where we accept responsibility under this clause, we shall at our sole option, repair or replace (as the case maybe) the items concerned which are proved to our satisfaction to have been lost or damaged prior to delivery to the customer.

## 8. Warranty

8.1 We warrant the goods or where we have erected or installed the goods, we warrant the installation and erection for a period of 12 months from the date of delivery or installation and erection ("the Warranty Period").

8.2 Our obligation under this warranty is limited, at our option, to repairing, replacing on

## 07. Index/Terms & Conditions · Register/AGBs

---

an exchange basis, any goods which are delivered with or develop defects in design, materials and workmanship under normal and proper use within the Warranty Period.

8.3 In the event of the customer becoming aware of the defect in the goods, or in the erection or installation during the Warranty Period, the customer shall promptly supply us within the Warranty Period with written particulars of such defect and use its best endeavours to provide us with all necessary access and other reasonable facilities and information and particulars required to enable us to inspect and remedy such defect.

8.4 In the case of damage in transit, the customer shall act in accordance with the provisions at clause 7.2 above.

8.5 If required by us, the customer shall return the defective goods to us, or at our option, in the case of goods which have been purchased by us from a supplier named in the quotation, to that supplier and

8.6 The customer shall bear the risk and cost of transport of the defective goods and of the repaired goods or goods supplied in replacement. If the complaint relating to the goods is justified we shall, at our discretion, reimburse the reasonable transport costs incurred by the customer.

8.7 We shall be under no obligation whatsoever to repair, replace, or make good any loss or damage or defect where such loss, damage or defect arises from the customer's neglect, misuse, faulty maintenance or installation or from alterations carried out without our prior written consent or from repairs carried out improperly by the customer or arising from normal wear and tear.

8.8 No setoff or deduction can be made from any account rendered by us for repairs. Where any damage has resulted from misuse or incorrect installation of the goods by the customer, we reserve the right to render an invoice for the amount of our charges based on our normal labour rates for the time spent and cost of materials or goods supplied.

8.9 Where any defective goods are replaced, the provisions of this clause shall apply to the replacement goods for the unexpired balance of the Warranty Period.

8.10 The customer shall not return any of the goods to us without our written consent and we shall not be under any liability whatsoever for the goods returned by the customer without such consent. If any goods are returned to us with our consent, we reserve the right to apply a handling charge. Any goods that are returned to us with our consent shall be accompanied by such certificates of decontamination as we may specify.

8.11 All goods supplied but not manufactured by us, are sold subject to the conditions of the sale of the manufacturer thereof and our sole liability in respect thereof shall be to give to the customer such benefits as we shall receive under any contract which we have with such manufacturer or under any guarantee which might be given to us in respect thereof. In the event of any failure by such manufacturer for whatever reason to accept such liability which may arise by reason of any defect in any product or part thereof we

shall be under no liability to the customer by reason thereof.

8.12 We shall not be liable for failing to perform the contract either wholly or in part if the failure is caused either wholly or in part by any circumstances outside our reasonable control.

## 9. Limitation of Damages

9.1 SAVE AS PROVIDED in the warranties set out in condition 8, we shall be under no liability to the customer for any damages or losses direct or indirect resulting from defects in design, materials, workmanship, installation or from any act or default by us whether negligent or otherwise.

9.2 We shall have no liability for any indirect or consequential losses or expenses suffered by the customer, however caused, including but not limited to loss or anticipated profits, goodwill, reputation, business receipts or contracts or losses or expenses resulting from third party claims.

9.3 Except in the case of death or personal injury caused by our negligence, or the negligence of our employees' agents or sub contractors, our liability under or in connection with these terms and conditions whether arising in contract, tort, negligence, breach of statutory duty or otherwise shall not exceed the contract price in respect of any occurrence or series of occurrences.

## 10. Property in the Goods

10.1 Notwithstanding delivery and the passing of risk in the goods or any other provision of these terms and conditions, the property in the goods shall not pass to the customer until we have received in cash, or cleared funds, in full all sums due from the customer to us under any agreement, transaction or series of transactions.

10.2 Until such time as the property and goods passes to the customer, the customer shall hold the goods as our fiduciary agent and bailee and shall keep the goods separate and apart from those of the customer and any third party and shall keep such goods properly stored, protected and insured and identified as our property.

10.3 Until such time as property and the goods passes to the customer, we shall have the right at any time to require the customer to deliver up the goods and, if the customer fails to do so forthwith, to enter upon any premises of the customer or any third party where the goods are stored with or without notice at any time to retake possession of the whole or any part of the goods to the value of all sums due to us without prejudice to any other right or remedy of ours.

## 11. Cancellation

No Order which has been accepted by us maybe cancelled by the customer except with our agreement in writing signed by a Director of ours and on terms that the customer

## 07. Index/Terms & Conditions · Register/AGBs

---

shall indemnify us in full against all loss (including loss of profits) (costs including the cost of all labour used, materials used) damages, charges and expenses incurred by us as a result of the cancellation.

### 12. Tools

Tools, moulds or other devices incorporating our know-how, which have been produced by us or handed over to the customer, remain our property, even if the customer has paid the costs for these in part or in full.

### 13. Statutory Limitation

All claims made by the customer - for whatever legal reasons - are subject to a limitation period of 12 months.

### 14. Validity of the Agreement

Should one or more of the provisions of this Agreement be or become invalid for any reason, the provisions should be interpreted to ensure that the commercial aim of the original, invalid provision is upheld. The validity of the remaining provisions will remain unaffected.

### **Edition 2015**

Poulten & Graf Limited  
Alfred's Way, Barking  
Essex. IG11 0AS, UK

## **Poulten & Graf GmbH – Conditions of Sale**

### 1. General conditions

Our General Terms of Business are binding for all deliveries and services. We will not accept general terms of business of our contract partners, which differ from our terms of sale, even without explicit objection. Additional agreements have to be confirmed by us in writing. Only German Law and German Text will be applicable. Place of performance is Wertheim. The place for jurisdiction, also for bill transactions, is Mosbach. All important data necessary for business transactions may be stored on computer. All information, which is not public knowledge and which results from the business relationship with our company, must be kept secret from third parties.

### 2. Offers

Our offers are always without engagement. In order to be valid, all transactions, orders and other agreements require our confirmation in writing. Our written order confirmation or agreement shall determine the nature and scope of the delivery. We reserve the right to make technical alterations to the offered products.

### 3. Risk, conditions of delivery and non-acceptance of delivery

The risk of ownership will be transferred to the buyer, as soon as the goods leave our company even in the event of delivery carriage-paid. Complaints can only be taken into account if submitted within 8 days of receipt of the goods. Delivered goods will only be accepted for return in their original packing. Illustrations and brochures etc. are not binding. Measurements and utilisation data are only to be regarded as approximate and are not binding.

Deliveries will be made ex factory, packing and insurance excluded, if no other agreement exists. Packing material of any kind may be returned to us within the range of legal regulations applicable to us. We will not pay for the cost of return of packing materials. In the event that dispatch or delivery is delayed at the request of the customer or for reasons for which he is responsible, the risk shall pass to the customer for the period of the delay, however at the request of the customer, we will arrange insurance cover that he requires at his expense.

Partial deliveries are permissible. Surplus quantities within the normal range (up to 15 %) are permissible and must be accepted and paid for by the customer. In case of delayed acceptance of the goods, we may sell the items at the customer's expense to a third party or we may ask for compensation.

### 4. Delivery times and Unforeseen Circumstances

Indicated delivery times are without obligation, unless a binding delivery time has been fixed in writing. Deliveries may be delayed due to Act of God or strikes, interruption of work, shortage of raw materials or failure to deliver by our suppliers. Liability for delayed delivery will be limited to proven, gross negligence, caused by us.

# 07. Index/Terms & Conditions · Register/AGBs

---

## 5. Special productions

Products, which are manufactured according to drawings, samples or instructions of the customer, cannot be returned. This includes in particular individual software for computers, delivered by us. All kinds of software delivered are subject to copyright and may be neither copied, modified nor presented in public without our prior written permission. Additional quantities within the normal range (up to 15 %) are permissible and must be accepted and paid for by the customer, unless exact quantities have been agreed upon in writing. The customer is responsible for ensuring that products manufactured according to his instructions, do not interfere with any trade right, copyright or other rights of a third party. The customer must compensate all damages resulting from any violation of such rights.

## 6. Payment conditions, prices, delayed payment

Products will be invoiced in EURO at the prices in effect on the day of delivery. Prices are ex factory excluding packaging and transportation insurance with immediate payment, unless otherwise provided. We reserve the right to include additional payment conditions for new customers. The customer may only set off payment against counter-claims, which are undisputed or legally determined. In case of delayed payment, we reserve the right to charge the customer with all costs as well as interest at 3 % above the respective EURIBOR 3 Month interest rate.

## 7. Warranty Claims

Warranty claims shall be recognised, if the defect has been reported to us in writing within 12 months of the delivery date. To the exclusion of further claims – subject to section 8 – we will provide a warranty as follows:

7.1 We shall have the option to either replace or to carry out repairs free of charge on all those parts which prove to be defective due to circumstances which occurred prior to the transfer of risk. The discovery of such defects must be reported to us immediately in writing. Parts replaced will become our property.

7.2 If we have to supply according to drawings, specifications, samples etc. provided by the customer, then the customer carries the risk of suitability for intended purpose. The deciding factor for the condition of the goods conforming to the agreement is the point in time of the transfer of risk according to section 3.

7.3 If a final inspection of the goods has been agreed or an inspection of the first sample, then customer's complaints will not be recognised later concerning defects which the customer should have noticed at the final inspection or inspection of the first sample if it had been carried out with due care and attention.

7.4 We must be given an opportunity to confirm the defect the complaint is being made about. Rejected goods must be returned to us immediately on demand. We will bear the transport costs if the customer complaint is justified. If the customer does not comply with these obligations or makes changes to the rejected goods without our agreement, then he forfeits any claims under warranty of quality.

7.5 We will bear the costs for replacement of goods including transport and reasonable costs for disassembly and assembly incurred due to the rectification of defects or the replacement delivery, providing the complaint proves to be justified.

7.6 Within the framework of legal provisions, the customer has the right to withdraw from the contract if, taking the legal exceptions into consideration, we do not successfully make use of a reasonable period of time given to us for rectification of defects or a replacement delivery due to a defect. If the defect is not significant, then the customer only has the right to claim a reduction of the contractual price. In all other cases the right to a reduction of the agreed price is excluded.

7.7 No warranty will be given in the following cases in particular: Unsuitable or inappropriate use, faulty installation or start-up by the customer or a third party, normal wear and tear, incorrect or negligent treatment, maintenance not carried out correctly, chemical, electrochemical or electrical influences providing we are not liable for them.

7.8 If the customer or a third party reworks the goods incorrectly, then we are not liable for the consequences arising from this. The same applies to changes made to the delivered item without our prior agreement. 7.9 We only provide a warranty for accessories added on within the framework of the warranty provided by our supplier.

## 8. Liability

8.1 If we are responsible for the customer not being able to use the delivery item as stipulated in the contract as a consequence of suggestions made and consultations carried out before or after the contract was concluded which were omitted or carried out incorrectly or through infringing upon other subsidiary obligations - in particular instructions for operating and maintaining the delivery item - then to the exclusion of further claims made by the customer, the stipulations in sections 7 and 8.2 apply correspondingly.

8.2 For damages, which have not been incurred by the delivered item itself, we only accept liability - whatever legal reasons are given - in the case of

- intent,
- gross negligence by executive employees,
- culpable injury of life, body, health,
- malicious silence with regard to a defect or where a defect's non-existence had been guaranteed
- defects of the delivery item for which liability must be accepted according to German Law on Product Liability for damages to persons or property with regards to privately used objects. In the case of culpable infringement of essential contractual obligations, we also accept liability for gross negligence of non-executive employees and in the case of slight negligence, in the latter case limited to reasonable, foreseeable damages typical to a contract.

Any further claims shall be excluded.

# 07. Index/Terms & Conditions · Register/AGBs

## 9. Reservation of proprietary rights

The delivered goods remain our property until complete payment. If the goods are resold or processed by the buyer, we are to be considered as manufacturer according to § 950 BGB and achieve the property of all intermediate or final products. The processing party is only the depositary.

For other items, which are not our property, we achieve the co-property of this new product in proportion to the goods reserved to the other components. The goods may only be sold within ordinary and regular business transactions and only, if claims arising from resale have not been assigned to a third party. The claims of the buyer resulting from resale are to be considered as assigned to us upon conclusion of the contract of sale. This is also valid, if our goods are combined or processed with other items. In this case, the assigned claims serve as our security only to the extent of the value of the merchandise sold under reservation of proprietary rights. The buyer is obliged to provide us with the names of third party debtors upon request and to inform them about the assignment. Amounts, which have been collected by him, have to be paid over to us, if our claims are due for payment. Distraints and hypothecations of the goods under reservation of proprietary rights or of the assigned claims are not permissible. The buyer has to inform us immediately about any access of a third party to the goods under reserve of proprietary rights or to the claims assigned to us. In the event of sales to foreign countries, should the retention of title not be admitted with the same effect as under German Law, the goods involved will remain our property until payment of all our claims resulting from the contractual relationship brought about by sale of said goods. In case, however, this retention of title is not admitted with the same effect as under German Law, but if it is allowed to reserve other rights to said goods, we will be entitled to exercise all of these rights. The buyer will be obliged to co-operate in all measures that we take for the protection of our right of ownership or, in its place, of any other title to the goods.

## 10. Tools

Tools, moulds or other devices incorporating our know-how, which have been produced by us or handed over to the customer, remain our property, even if the customer has paid the costs for these in part or in full.

## 11. Statutory Limitation

All claims made by the customer - for whatever legal reasons - are subject to a limitation period of 12 months. For wilful or fraudulent behaviour and for claims based on the German Product Liability Law, the legal limitations apply.

## 12. Validity of the Agreement

Should one or more of the provisions of this Agreement be or become invalid for any reason, the provisions should be interpreted to ensure that the commercial aim of the original, invalid provision is upheld. The validity of the remaining provisions will remain unaffected.

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